



Medical Diode Laser Hair Removal

1800W+1200W handles


Why all the beauty salons & beauticians choice titanium

1. Professional. we have more than 14 years experiences in the High-end beauty machine area.
2. Safety. no side effect on skin and sweat glands, no scar.
3. TEC refrigeration system, fast cooling, which achieve painless hair removal.
4. Popular. Triple wavelength: 755nm+ 808nm+1064nm ,suitable for all the skin types, black skin fine hair all can be easily removed.
5. High energy. Lager spot size, 20-40 million shots, effectively destroy hair follicles, permanent hair removal. Also prompt treatment speed and efficiency.

Handle show






Laser show Handle laser bar optional



⊙ 10bar


⊙ 12bar

⊙ 16bar



Working Principle

Suits for all skin and hair types




- **755nm for white skin** (fine, golden hair removal)
- **808nm for yellow neutral skin** (brown hair removal)
- **1064nm for black skin** (black hair removal)



Customer Feedback

Treatment Comparison



**Please consult us for more details.
24 hours online!**