



Descrizione dei prodotti

Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



Laser a diodiElight/IPLNd yagRF, 4 diverse impugnature in un'unica macchina

1. Questo laser 4 in 1 RF Elight IPL laser a diodi nd yag ha accessori di alta qualità, forte potenza *
La maniglia del laser a diodi utilizza la barra laser dorata importata dalla Germania, 50 milioni di scatti.


*La maniglia Elight utilizza la lampada del Regno Unito, 150.000-30.000 scatti

*Nd yag: equipaggia 5 testine di trattamento, regolabili 1064 nm/532 nm, fisse 1064 nm/532 nm, fisse 1320 nm.


*RF: equipaggia 6 testine di trattamento, utilizzate per la rimozione delle rughe, il lifting della pelle, il rassodamento della pelle.

Three wavelengths work together

12*24mm² {




12*35mm² {



Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.





ELIGHT HANDLE

specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15*50mm, 15*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm
Standard: 420/510/680nm



Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

HANDLE-ADVANTAGE



specification

- **Multiple choice**
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- **Multifunction**
1320nm for skin rejuvenation
- **Indicator light**
Dot position is the indicator light position dot more accurate



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm



Fixed
1064nm



Before and after



RF Handl e



12mm
eye-rf

28mm small
tripolar rf head

36mm big
tripolar rf head

14mm
bipolar-rf

28mm
bipolar-rf

38mm
bipolar-rf

BEFORE & AFTER ^R



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

Advantage



808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow, embroidery removal, freckle removal.

RF handle

Skin tightening/lifting, wrinkle removal on body

**Fo ulteriori informazioni sull'introduzione del prodotto e sui dettagli sull'uso, consultateci immediatamente
Online 24 ore!!**