



# **MULTI-FUNCTION 4 IN 1 MACHINE**

Descrizione dei prodotti

# Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



## **Laser a diodi Elight/IPL Nd yag RF, 4 diversi manici in un'unica macchina**

1. Questo 4 in 1 RF Elight IPL laser a diodi laser nd yag ha accessori di alta qualità, forte potenza \*  
L'impugnatura del laser a diodi utilizza la barra laser dorata importata dalla Germania, 50 milioni di scatti.


\*La maniglia Elight utilizza una lampada del Regno Unito, 150.000-30.000 scatti

\*Nd yag: Equipaggia 5 testine di trattamento, regolabili 1064nm/532nm, fisse 1064nm/532nm, fisse 1320nm.


\*RF: Equipaggia 6 testine di trattamento, utilizzate per la rimozione delle rughe, il sollevamento della pelle, il rassodamento della pelle.

## *Three wavelengths work together*

**12\*24mm<sup>2</sup>** {




**12\*35mm<sup>2</sup>** {



**Suitable for hair removal on any part of the body,  
Such as arms, legs, face, armpits, back, chest, etc.**





# ELIGHT HANDLE

## specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15\*50mm, 15\*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm  
Standard: 420/510/680nm



## Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

# HANDLE-ADVANTAGE



## specification

- Multiple choice  
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction  
1320nm for skin rejuvenation
- Indicator light  
Dot position is the indicator light position dot more accurate



Adjustable  
532nm



Adjustable  
1064nm



Black doll  
1320nm



Fixed  
532nm



Fixed  
1064nm



# Before and after



## RF Handl e

12mm eye-rf	28mm small tripolar rf head	36mm big tripolar rf head	14mm bipolar-rf	28mm bipolar-rf	38mm bipolar-rf

# BEFORE & AFTER



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

## Advantage



### 808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

### Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

### ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow, embroidery removal, freckle removal.

### RF handle

Skin tightening/lifting, wrinkle removal on body

**Fo ulteriori informazioni sull'introduzione del prodotto e sui dettagli di utilizzo, si prega di contattarci immediatamente  
In linea 24 ore!!**