

EXFU Seiko



11D HIFU

Macro & Micro High Intensity Focused Ultrasound System



SAFE



EFFECTIVE



MULTI-DEPTH



EXFU Seiko 11D HIFU is a non-invasive, non-surgical treatment for skin rejuvenation. It uses high-intensity focused ultrasound energy to stimulate collagen production and tighten the skin. The device is safe, effective, and multi-depth, allowing for treatment of different skin layers. It is suitable for various skin types and conditions, including wrinkles, sagging skin, and loss of elasticity. The treatment is performed in a clinical setting and typically takes 30-60 minutes. The results are visible immediately and continue to improve over time. For more information, please contact your healthcare provider.





Macro & Micro High Intensity Focused Ultrasound System

7CARTRIDGES+4 GUN HEADS

1.5mm,2.0mm,3.0mm,4.5mm for face treatment
 6.0mm.9mm.13mm for body treatment.

此系统采用先进的超声技术，能够实现精准的能量聚焦，对皮肤深层组织进行非侵入性治疗。系统包含7种不同规格的超声探头，适用于面部及身体不同部位的紧致和塑形。其独特的能量分布模式能够有效刺激胶原蛋白再生，改善皮肤弹性和紧致度。系统操作简便，治疗过程舒适，无恢复期，是您进行面部轮廓提升和身体塑形的首选方案。







MF2.0MM

Unique ultra-thin
probe Clexis patented
technology

"Special cartridge for eyes treatment"



Body Treatment



• Treatment Effect •



□□□□□	7 □□□□ □□□□ □□□□□□□□□□
□□□□	□□□□ □□□□ □□□□ □□□□□□□□□□

□□□□□	□"□ 4.5 ,□"□ 3.0 ,□"□ 1.5 :□□□□ □□□□
□□□□□ □□□□□	7D HF: MF2/6/9/13 L4-4.5 L7-1.5 L7-3.0
□□□□□ □□□□□	0.1J-2.5J (0.1J □□□ □□□)
□□□□□ □□ □□□□	1.0 □"□ 4.0-□"□
□□□□□ □□□□	5 □"□ 25-□"□
□□□□□	□□□□□ 300000 □□□□□ □□
□□□ □□□□	15 □□□□'
□□□□	110V-240V, 50/60Hz



- □□ □□ □□ □□□□□
- □□□□□ , '□□□□ 15.6 □□□□□ □□□□ □□ □□□ .1
- .□□□□□□□ □□ □□□□ 6 ,□□□ □□□□ ,□□□□/□□□ .2
- .□□□ □□ □□□□□□□□□□ .3
- □□□□□ .4 OEM □□□□□/□□□□□) □□□□□...□
- □□ □□□□□□ □□□□□ □□ □□□□□□□ □□□□□□!