



GENES LASER

NEWEST COMING, FASTEST
STRONGEST, DOUBLE HAIR
REMOVAL SYSTEM

1200W + 1800W TITANIUM COMING

VIEW MORE



- * 1800 000000w 000000 0000 000000 ,000000 000000 000000000000 00 000000 0000000!
- * 00000 0000000 .0000 00000000 00000 000000 00000 00000 00000 00000 00000 00000 00000 000000 000000 00000 00000000 000000 0000 ,00000 00 000000000000 00000000 00000000 00000000 00000 00000 6 0000 00000000 00 00.
- *00 ...00000000 ,00000000 ,00 ,000 ,000 0000 ,00000000 ,00000 0000 0000000000 00000 00 00000 0000 000000 00000 000000 00 00.
- *1-10 000hz.0000000 0000 .00000000 0000000 00000 0000000 000000 !!!00000 000000!!

Handle display

① Optional 12*24mm/12*36mm light spot size.

Different spot size meet different body parts hair removal needs .



spot size : 12*24mm



spot size : 12*36mm

② 10bars 1200w / 15bars 1800w. 10million—40million shots.

Effective for permanent hair removal

③ 755nm and 808nm and 1064nm , suits for all skin and hair types.

Cooling effect

The handle has a cooling effect that improve treatment comfort.



spot size : 12*24mm



spot size : 12*36mm



755 nm 808nm 1064nm, suits for all skin and hair types


755nm .(epidermal/dermal) hair removal

808nm (epidermal) hair removal

1064nm .() () () ()
() () () () ()

working principle


-  755 nm for white skin (fine, golden hair removal);
-  808nm for yellow neutral skin (brown hair removal);
-  1064nm for black skin (black hair removal)



A natural




B energy absorbed by hair



C hair follicle died

working principle


Triple wavelength laser spots work together



□□□□ □□□□ □□□□ □□□□□□

Scope of application

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...



CE and TUV CERTIFICATE

**U.S. FOOD & DRUG
ADMINISTRATION**

Home | Food | Drugs | Medical Devices | Radiation-Emitting Products | Vaccines, Biologics

510(k) Premarket Notification

FDA Home | Medical Devices | Databases

510(k) | Definitions | Registration & Listing | Adverse Events | Recall | CFR Title 21 | Radiation-Emitting Products | X-Ray Assemblies

Search

Device Classification Name	External Laser Box
510(k) Number	K200525
Device Name	Medical Grade Laser
Applicant	Welfang Mingliang & Co., Ltd., Yuyang 5, Wufang, CN 261000
Applicant Contact	Li Mingliang
Correspondent	Beijing Believe Med. Co., Ltd., Building 1, Fangshan District, Beijing, CN 102401
Correspondent Contact	Ray Wang
Regulation Number	823.6533
Classification Product Code	GE2
Date Received	03/02/2009
Decision Date	05/27/2009
Decision	Substantially Equivalent
Regulation Medical Specialty	General & Plastic Surgery
510k Review Panel	General & Plastic Surgery
Type	Traditional
Reviewed By Third Party	No
Combination Product	No

Certificate of Compliance

No. 04170333-WFSCY02
Technical Construction No. 10, 14170333-WFSCY02

Certificate's Holder: Welfang Mingliang Electronics Co., Ltd.
Floor 3, Workshop 2, No.13426, Yuyang 5th Street, Yuyang Community, Xincheng Street, High-tech Zone, Wufang, Shandong, China.

Certification ECM Mark: 

Product: LASER BEAUTY DEVICE
Model(s): (See the following annex)

Verification to: Standard: EN 60601-1-2:2012 EN 60601-1:2006+A2:2014
Issued to CE Directive(s): Medical Devices (Directive 93/42/EEC, 2007/45/EC)

Remarks: The products have been verified on a technical basis. The products comply with the requirements of the CE Directive(s) in relation to the above listed standards. The above Compliance Mark can be affixed on the products according to the CE regulation about its position and its use. The regulation can be found at www.ecmcert.com. The Certificate of Compliance can be checked for validity at www.ecmcert.com.

Additional Information: Application issued to CE mark the product EN 60601-1-2:2012 EN 60601-1:2006+A2:2014. The manufacturer is responsible to ensure the CE marking Certificate is valid through an approved notified body and the products of the following activities are marked by the CE mark as required by the notified body. Before putting the CE mark on the products.

CE

Date of issue: 01 December 2013
Date of expiry: 31 December 2015

ECM
11, Welfang Mingliang Electronics Co., Ltd. (www.welfangmingliang.com.cn)
www.ecmcert.com

ITC
Management Systems Certification Body No. 2022
Institute for Testing and Certification, Inc.
1, Bell 208, 754 21 Zhen, Czech Republic
www.itcinc.com

CERTIFICATE

No. 10 0144 SJ

We certify on the basis of a performed audit that company
Welfang Mingliang Electronics Co., Ltd.
338Wong Fong Dongxue Street Wufang China
Company Reg. No.: 3707602000018

has implemented and documented a functional quality management system
in compliance with the requirements of the standard
EN ISO 13485:2003

Covering the following activities:
Design & Development, Production and Distribution of IPN, Phototherapy System,
MD, YAG Q Switched Laser and RF Device

The Certificate is issued on the basis of the results mentioned in Audit Report No.
20090802017. The Certificate validity is conditioned to positive results of surveillance audits
which the certified company conducts at regular intervals.


Date of issue: 01.12.2009
Valid until: 31.12.2010
Date of the last certification audit: 10.07.2007

Paul Vlas
Eng. Pavel Vlas
Head of Certification Body



SHIPPING

Welcome to visit our factory



In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month