



# Medical Diode Laser Hair Removal

1800W+1200W handles

1. 15-20 minutes per session, 1-2 sessions per week.

2. 1-2 sessions per week, 1-2 sessions per week.

3. 1-2 sessions per week, 1-2 sessions per week.

4. 755: 1-2 sessions per week, 1-2 sessions per week.

5. 1-2 sessions per week, 1-2 sessions per week, 50, 1-2 sessions per week.

1-2 sessions per week, 1-2 sessions per week, 50, 1-2 sessions per week.

1-2 sessions per week

# 12\*24mm or 12\*36mm handle

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...



- >>> Big 20\*12 □□□□□ □□□□mm, 12\*35mm for option
- >>> Laser technology, big power 10bar 1200W, 12 bar 1800W
- >>> 300g handle, light weight, free quick slide treatment



>>> New nose hair treatment head, applicable to small part treatment  
>>> 1 minute refrigeration frosting, painless treatment experience






ALEX 755nm  
**For white skin**  
(fine, golden hair removal )



SPEED 810nm  
**For yellow neutral skin**  
(brown hair removal)




YAG 1064nm  
**For black skin**  
(black hair removal)



## Treatment Effect

This is the clinical result of hair removal  
It takes 3--5 times to realize  
permanently hair removal





- \* 4 ק"מ 14 סממנים
- \* OEM סממנים
- \* סממנים

### מפרט של מכונה להסרת שיער בלייזר דיודה

פריט	755nm/808nm/1064nm לייזר דיודה אנכי 755
סוג לייזר	בר לייזר רפואי דיודה לייזר: 10 בר/ים/15 בר/ים
כוחו של כל בר	120w/Bar רפואי עם
מערכת תצוגה	'מסך מגע בצבע אמיתי בגודל 14 אינץ'
אורך גל	808nm/755nm/1064nm
Spot Size	24*12mm <sup>2</sup> /12*36mm <sup>2</sup> עם רפואי
תוחלת חיים של ידית לייזר	לפחות 40 מיליון יריות
אנרגיה	(אופציונלי J/cm <sup>2</sup> /מתכוון) עד 150 J/cm <sup>2</sup> 10-150
משך הדופק	מתכוון 30 ~ 400ms
תדירות	1.0-10 הרץ;
מערכת קירור	קומפרסור אוויר סגור זרימת מים
ידיית כוח	1000w/1200w/1600w/1800w
כוח מכונה	וואט 4000
מקור חשמל	220V/50HZ 110V/60HZ או
טמפרטורת קריסטל	מעלות צלזיוס 3~20-
טמפרטורת המכונה	מעלות צלזיוס 24
דרישת חשמל	AC220V/110V 50~60Hz

□□□□□□ □□□□

# Free Shipping

## Custom Logo

## Hot Selling Big Spot





# Medical diode laser Hair removal machine

V19



• 100% 100% 100% 100% 100%  
100% 100% 100% 100%  
100% 100% 100% 100% 100% 100%  
100% 100% 100% 100%  
100% 100% 100% 100% 100%  
• 100% 100% 100% 100% 100%  
100% 100% 100% 24!