



□□□□□ □□□□□

Real shot display



808 handle



Elight handle



ND YAG Laser handle




RF handle




1-4 RF Elight IPL * 50,000-30,000, 1064nm/532nm, 1320nm.
 *Nd yag: 1064nm, 532nm, 1320nm.
 *RF: 1-4 RF Elight IPL * 50,000-30,000, 1064nm/532nm, 1320nm.

Three wavelengths work together

12*24mm² {




12*35mm² {



**Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.**





ELIGHT HANDLE

specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15*50mm, 15*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm
Standard: 420/510/680nm



Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

HANDLE-ADVANTAGE



specification

- Multiple choice
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction
1320nm for skin rejuvenation
- Indicator light
Dot position is the indicator light position dot more accurate



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm



Fixed
1064nm



Before and after



RF Handl e

12mm eye-rf	28mm small tripolar rf head	36mm big tripolar rf head	14mm bipolar-rf	28mm bipolar-rf	38mm bipolar-rf

BEFORE & AFTER [®]



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

Advantage [®]



808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow embroidery removal, freckle removal.

RF handle

Skin tightening/lifting, wrinkle removal on body

For more information, please contact us at 24 hours!!