

**DIODE LASER**



**ICELLO**



[illegible]

50!




120W/ bar 1800W  
Spot size :12\*36mm



120W/ bar 1200W  
Spot size :12\*24mm


TREATMENT HEAD  
REMOVING NOSE HAIR



Eye

Ear

Nose



Mouth

Hand

Diameter 8mm



// Hair removal on multiple parts



1. 32-°C 温度点 温度点 温度点 .温度点 温度点°C°, 温度点 温度点 温度点 温度点
2. 温度点 温度点 温度点 温度点 温度点 温度点 2-3 .温度点 温度点 温度点
3. 温度点 温度点 温度点 温度点 温度点 温度点 温度点 温度点 温度点
4. 温度点 温度点 ,温度点 温度点 温度点 ,温度点 温度点 温度点 温度点 温度点 温度点 温度点 温度点 .温度点 温度点 温度点 温度点 温度点.



755/808/1064 nm 温度点 温度点 温度点 温度点 nm.

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>America Imported Emitter, strong energy

>>10/12 pieces or 16 pieces Germany Dilas original laser bars,1200W/1800W power is capable to meet all demands of hair removal on the whole body


# Medical diode laser

**50000000+shots**  
**can treat 10000+clients**

- \*Not easy to be burned
- \*Low maintenance cost
- \*Stable quality



- >>1-10Hz fast frequency,faster, smooth and effective treatment
- >>Super sapphire touch cooling system ensures painless and comfortable treatment
- >> Italy imported stainless water pump, it can keep machine working long time without stop. Make sure the heat were removed very quickly
- >>Double filters to keep water quality





## The real experience is visible



Through live test verification ,safe ,long,painless  
and effective



Before



After



EXFU®  
Simple · Safe · Enjoy

**14 inch**  
**4K HD screen**  
Ultra thin technology version



4 '0000 14 000K, 0000 000 00000

\* 6 000000 000 0000000

00000 000000/00000 0000000000 000 000000 00 \*

0000 00 000000 0000 0000 0000000 00000 \*

## Customized Service

● Grey



● Purple



● Black



● Gold



# You'll get

Two major service guarantees  
Worry-free after-sales

## 01 Professional after-sales group

24 hours offers the solution  
Solve problems quickly with your active cooperation

## 02 Free maintenance

Warranty period, Send your replacement  
parts for free

Our products have passed 15 years of strict testing and are certified FDA, CE, ISO13485.  
We guarantee the quality of our products and provide you with the best service!

## Establishment Registration & Device Listing

FDA Home Medical Devices Databases



New Search

Back To Search Results

### Establishment:

WEIFANG MINGLIANG ELECTRONICS CO.,LTD.  
Floor 3, Building 4, Second Accelerator, Yuqing Eastern Street, High-Tech District  
Weifang Shandong, CN 261000  
**Registration Number:** 3007549439  
**FEI Number\*:** 3007549439  
**Status:** Active  
**Date Of Registration Status:** 2022

### Owner/Operator:

Weifang Mingliang Electronics CO.,LTD.  
Floor 3, Building 4, Second Accelerator, Yuqing Eastern Street, High-Tech District  
Weifang, Shandong CN 261000  
**Owner/Operator Number:** 10082104

### Official Correspondent:

Pingping Lv  
Floor 3, Building 4, Second Accelerator, Yuqing Eastern Street, High-Tech District  
Weifang, Shandong CN 261000  
**Phone:** 86-10-56335780

### US Agent:

Wen Xu  
4611 Maydew Mews Way Ellicott City MD  
Ellicott, MD US 21042  
**Phone:** 443 8968583 Ext  
**Email:** Infor@Believe-Med.Com

\* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English



### Authentication

FDA,NQA13485,CB certification  
and more Complete certification  
materials



### 60000+ Customers

We have professional after sale  
team and we provide 7\*24  
services



### 2 Years Warranty

2 years warranty for each  
machine, All questions will get  
reponse within 24 hours

7 00000 00000 !0000000 0000000X24 0000!