




Soprano ice PLATINUM

Medical beauty machine

Medical diode laser hair removal
For all skin types,
triple wavelength laser work together



.000000 000000 0000 000000 000000
;000 00 000 000 000000 .0000000000 000*
;000 00 000 000 000000 *
.000000 0000 *
;0000 00000000 0000000000 0000 0000 *
808 0000000000 00000000 00 0000 *nm 7558081064 00nm;

350g

Ultra light handle



12×36mm

Handle

1800W diode laser hair removal



12×24 mm


Handle

1200W diode laser hair removal

00000 **Three Waves 755** 0000nm 808nm 1064nm, 0000000 0000 0000 000 00000000 000000.

755nm .(00000/00 0000) 000 0000 0000 000000

808nm .(목발) 목발 목발 목발 목발
 1064nm .(목발 목발) 목발 목발 목발 목발
 목발 목발 목발 목발 목발 목발 목발 목발 목발 목발




ALEX 755nm
For white skin
 (fine, golden hair removal)



SPEED 810nm
For yellow neutral skin
 (brown hair removal)



YAG 1064nm
For black skin
 (black hair removal)



Treatment Comparison

