

Medical



1200W+1800W high power handle
14-inch 4K HD large screen
755nm+808nm+1064nm
Painless permanent hair removal

Diode Laser

□□□□□□ □□□□□□ □□□□□□ □□□□ □□□□□□ □□□□ □□□□□□

4 '□□□□ 14 □□□□□□ □□□□ □□□□K HD

1200 □□□□□□ □□W1800W □□□□□□ □□□□□□

□□□□ □□□□□□□□□□ □□□ □□□□□□ □□□□□□

□□□□ □□□□ □□□□ □□□□ ,□□□□ 350-□ □□□□□ □□□□□□

□□□□□□ □□□□□□ □□ □□□□□□ □□□□ □□ □□□□□□ □□□□ ,□"□ 36 * 12 □"□ 24 * 12 □□□□□□

Machine Display



1. .0000 00000 00000 00000 0000 0000 0000 0000 ,0000 00000 ,0000 0000 0000
755 00000 000000 0000 0000 000000 .2nm, 808nm, 1064nm
3 00000 00000.TEC, 000 00 0000 000000



Before and After

Treatment comparison



□□□□□□ □□ ,□□□□□□ □□□□□ □□□□□□ .□□□□ □□□□ ,□□□□□□ ,□□□□ □□□□□ □□□□□ □□□□□ 100% □□□□□□!