



□□□□□ □□ □□□□□ □□□ .□□□□□□ □□□□□ □□□□□ □□ □□□□□ □□□ □□□□□ □□□□□ 24 □□□□□ □□□□□ !!

Diode Laser

For You. For Life.



1064 808 755

- 1) TEC, 24
- 40,000,000-2
- 1600

350g lightweight **handle**, more ergonomic

Available in **black or white**



755 808 1064 1000 1200 1600 20*12 35*12 35*12

1000 10,0" 20*12 1000w.

1200 12,0" 35*12 1200w.

1600 16,0" 35*12 1600w.

1000 1200 1600 20*12 35*12 35*12

1000 1200 1600 20*12 35*12 35*12 1000 1200 1600 20*12 35*12 35*12



12*24mm 1200W handle
Additional therapeutic head



12*36mm 1800W handle
Additional therapeutic head






755 808 1064 000000 000000 0000 000000 0000 0000

- 1) 000000 0000 000000 00000000 0000 0000
- 2) 0000 0000 000 000000 0000000000 0000 0000
- 000000 000 00000 000000 ,0000 00000000 0000 00000000 (3.



00000 0000 000 000 ,0000000 000000 00000








1. 14.1 英寸触摸屏
 2. 6 英寸触摸屏, 触摸屏/触摸屏
 3. 触摸屏
 4. OEM (触摸屏/触摸屏) 触摸屏...

The new **14-inch 4K** screen

Support system customization



产品特点

波长	755nm 808nm 1064nm
脉冲宽度	16 --- 8ms, 可调
光斑大小	100mm
治疗面积	14英寸触摸屏
冷却系统	808nm 808nm 755nm 1064nm
治疗头	12*35mm 20*12mm
治疗深度	808nm 10,000,000-40,000,000 脉冲
治疗能量	10-150 J/cm² 150 J/cm² (可调)
治疗时间	40 ~ 800ms
治疗次数	1.0-10 次
治疗模式	TEC 冷却系统
功率	800w/1000w/1200w/1600w
电压	2000W
频率	220V/50HZ 110V/60HZ
温度范围	-10~3℃
重量	27~31kg
电源	AC220V/110V 50~60Hz

产品优势, 产品优势!