



808 Diode Laser

Big Power: 1000w 1200w 1600w

Triple wavelength: 755nm 808nm 1064nm

For different skin type: white, brown ,black

□□□□□□ □□□□□□

Machine Parts



□□□□ □□
□□□□□□ □□□□

□□□□□□ □□□□□□ □□□□ □□□□□□ □□□□□□ □□□□□□
755nm/808nm/1064nm

Handle display Handle color optional

weight **450g**

The handle weighs only 450g, making the treatment process easier



Cooling effect

Fast cooling, more comfortable treatment process



12*20mm



12*35mm




12*20mm



12*35mm

TEC 24 24

Laser show Handle laser bar optional



1064nm: 1064nm 808nm 755nm
 755nm: 755nm 808nm 1064nm
 ND YAG 1064nm: 1064nm 808nm 755nm

Before and After

Treatment comparison



Product Display



!!! 00000000 ,000000 0000 000000!