




- 1.3 waves 1064nm + Diode 808nm + 755nm, 3 in 1 all skin type can be used
High power: NEW 1000W/ 1200W /1600W handle for options
2. Big spot size 12*20mm, 12*35mm face & body hair removal solution!
3. Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
4. Long life: double life support guarantee

Working Principle

Suits for all skin and hair types



- **755nm for white skin** (fine, golden hair removal)
- **808nm for yellow neutral skin** (brown hair removal)
- **1064nm for black skin** (black hair removal)





Different spot size Handle color optional



Safe Use On All Skin Type

Alex 755nm for light skin

Speed Diode 808nm for yellow skin

Long Pulse 1064nm for darker skin

Treatment effect



Software & Program

1. Male & femal, 6 body parts, 6 skin types setting smart software
2. 6 language options, your native language can be added
3. OEM program service available(logo/system...)

