

## **2021 Salon Queen**

**Trending Model 1600W Laser  
4K Smart Screen**

Newest Germany Tech  
Super Ice Cooling, Pain-free  
Get most comfortable treatment



Already help over customers from 102 countries make good business.

Various templates for you to choose.

Product Description




## **Advantages**


1. 3 waves Yag 1064nm + Diode 808nm + Alex 755nm, 3 in 1 all skin type can be used  
High power: NEW 1000W/ 1200W /1600W handle for options
2. Big spot size 12\*20mm, 12\*35mm, face & body hair removal solution!
3. Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
4. Long life: double life support guarantee



**A** natural



**B** energy absorbed by hair



**C** hair follicle died

### **Safe Use On All Skin Type**

Alex 755nm for light skin

Speed Diode 808nm for yellow skin

Long Pulse 1064nm for darker skin


# Handle Details

Three Wavelengths Work Togher




**12\*20nm**

755nm Laser  
808nm Laser  
1064nm Laser



**12\*35nm**

755nm Laser  
808nm Laser  
1064nm Laser



## Features & Diode Laser

1. Laser bar: Imported laser bar US Coherent/ Dilas brand
2. 3 waves: 755nm 808nm 1064nm 3 in 1
3. Spot size: 12\*20mm, 12\*35mm
4. Power: High power 1000W 1200W 1600W & OEM
3. Cooling: Super cooling start frozen in 1 minutes, 24 hr continue work



### **Pain-free Guarantee**

Super TEC cooling system Fast cooling 1 minute frozen

Lowest temperature -22°C



# Before and After

## Treatment comparison



### Software & Program

1. Easy Start: Based on male & female/6 body parts, Add 6 skin types setting smart software
2. 6 language options, your native language can be added
3. OEM program service available (logo/system...)

