

DIODE LASER



ICELLO



MLKJ is the best commercial laser hair removal machine.

According to the principle of selective photothermodynamics, the wavelength, energy and pulse width of the laser can be properly adjusted so that the laser passes from the skin surface to the root hair follicle. The light energy is absorbed and converted into heat energy that destroys the hair follicle tissue, making the hair unable to regenerate, making MLKJ lasers painless and permanent with patented technology. The laser hair removal machine is a non-invasive treatment method, after the treatment it will not cause any special damage to the surrounding skin, it will not affect the normal life and work.

□ FDA and medical CE certified device

□ 3 in 1 - three laser technologies in one device

□ Alexandrite, diode and Nd:YAG lasers - treat all types of skin

□ Easy operation

□ Online training is available

□ 50 million shots per handpiece!

Handle display



1. Without pain. The tip of the handpiece cools down to -32 degrees Celsius and makes the treatment comfortable and effective.
2. Primary results of treatment. 2-3 treatment sessions can show permanent results
3. Virtually painless in contact with sapphire, unparalleled comfort and safety
4. The tip of the face increases accessibility. Easily treat areas that are usually hard to reach, including ears, nostrils and glabella.



The price of laser hair removal machine 755/808/1064nm.

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>America Imported Emitter, strong energy

>>10/12 pieces or 16 pieces Germany Dilas original laser bars, 1200W/1800W power is capable to meet all demands of hair removal on the whole body


Medical diode laser

50000000+shots
can treat 10000+clients

- *Not easy to be burned
- *Low maintenance cost
- *Stable quality



- >>1-10Hz fast frequency,faster, smooth and effective treatment
- >>Super sapphire touch cooling system ensures painless and comfortable treatment
- >> Italy imported stainless water pump, it can keep machine working long time without stop. Make sure the heat were removed very quickly
- >>Double filters to keep water quality



The real experience is visible

Through live test verification ,safe ,long,painless
and effective



14 inch 4K screen, smart touch experience

- * 6 skin styles to choose from
- * Easy to provide operational skills / training services
- * Free logo and language can be customized on the screen



You'll get

Two major service guarantees
Worry-free after-sales

01 Professional after-sales group

24 hours offers the solution
Solve problems quickly with your active cooperation

02 Free maintenance

Warranty period, Send your replacement
parts for free

MLKJ

Manufacturer of beauty tools for 15 years with FDA, CE, ISO13485 certification.
We are committed to making your beauty business easier!

Establishment Registration & Device Listing

FDA Home Medical Devices Databases



New Search	Back To Search Results
<p>Establishment: WEIFANG MINGLIANG ELECTRONICS CO.,LTD. Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District Weifang Shandong, CN 261000 Registration Number: 3007549439 FEI Number*: 3007549439 Status: Active Date Of Registration Status: 2022</p> <p>Owner/Operator: Weifang Mingliang Electronics CO.,LTD. Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District Weifang, Shandong CN 261000 Owner/Operator Number: 10082104</p> <p>Official Correspondent: Pingping Lv Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District Weifang, Shandong CN 261000 Phone: 86-10-56335780</p> <p>US Agent: Wen Xu 4611 Maydew Mews Way Ellicott City MD Ellicott, MD US 21042 Phone: 443 8968583 Ext Email: Infor@Believe-Med.Com</p>	
<p>* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set</p>	

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English



Authentication

FDA,NQA13485,CB certification
and more Complete certification
materials



60000+ Customers

We have professional after sale
team and we provide 7*24
services



2 Years Warranty

2 years warranty for each
machine, All questions will get
reponse within 24 hours

Welcome to the consultation! 7x24 hours online service!