



Products Description

# Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



## **Diode laser+Elight/IPL+Nd yag+RF, 4 different handles into one machine**

1. This 4 in 1 RF+Elight IPL+diode laser+nd yag laser have high quality accessories, strong power\*Diode laser handle use Germany imported golden laser bar, 50 million shots time.


\*Elight handle use UK lamp, 150,000-30,000 shots

\*Nd yag :Equip 5 treatment heads, adjustable 1064nm/532nm, fixed 1064nm/532nm, fixed 1320nm.


\*RF :Equip 6 treatment heads, used for wrinkle removal, skin lifting, skin tightening.

## *Three wavelengths work together*

**12\*24mm<sup>2</sup>** {




**12\*35mm<sup>2</sup>** {



Suitable for hair removal on any part of the body,  
Such as arms, legs, face, armpits, back, chest, etc.







## ELIGHT HANDLE

### specification

- **High-quality:** Elight handpiece with super cooling system
- **Optional Spot Size** 15\*50mm, 15\*35mm
- **Long shots time:** 150000 times
- **7 filters:** 420/510/560/585/640/690/755nm  
Standard: 420/510/680nm



## Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

# HANDLE-ADVANTAGE



## specification

- **Multiple choice**  
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- **Multifunction**  
1320nm for skin rejuvenation
- **Indicator light**  
Dot position is the indicator light position dot more accurate



Adjustable  
532nm



Adjustable  
1064nm



Black doll  
1320nm



Fixed  
532nm



Fixed  
1064nm





## Before and after



## RF Handl e



12mm  
eye-rf

28mm small  
tripolar rf head

36mm big  
tripolar rf head

14mm  
bipolar-rf

28mm  
bipolar-rf

38mm  
bipolar-rf

# BEFORE & AFTER <sup>R</sup>



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

## Advantage



### 808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

### Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

### ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow, embroidery removal, freckle removal.

### RF handle

Skin tightening/lifting, wrinkle removal on body

**For more information about product introduction and use details, please consult us immediately  
Online 24 hours!!**