



1. Rotable handle plug, easier to install.

2. Wear-resistant tube, more durable. Handle is injection molding design, lighter and more practical.

3. Brand new 14 inch screen is clear and easy to operate.

4. Handle is injection molding design, lighter and more practical. 5. Upgrade the heat dissipation hole, the heat dissipation effect is good, and the heat dissipation is faster!



>>Germany Emitter, strong energy

>>1200W power of hair removal on the whole body;

>>1-10Hz fast frequency

>> Italy water pump, it can keep machine working 24 hours without stop.



Best cooling system:air +water+TEC semiconductor condenser□

Water temperature can be controlled between 26°C—28°C, 24 hours continuously working.

Working Principle

755nm and 808nm and 1064nm , suits for all skin and hair types.



ALEX 755nm

For white skin

(fine, golden hair removal)



SPEED 810nm

For yellow neutral skin

(brown hair removal)



YAG 1064nm

For black skin

(black hair removal)

Handle Long life time 20-40 millions shots

Speed 808nm: Golden standard wavelength for all hair types

Alex 755nm: Specific effective for blonde hair, fine hair

Long Pulse 1064nm: Specific effective for dark, tanned skin best quality opt laser/ nd yag /alexandrite laser hair removal machine price

Machine Display



Packing

safe and thoughtful packaging service



Built-in prevention Collision bubble

Two-year warranty for machine handle as consumable,
one-year warranty.

OEM&ODM service for distributor.

Tight packaging

Training: (Video+Manual+Online Service)
Can Keep You Operate It Easily.



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service,