



1. Three wavelength diode laser.
2. Golden laser bars from German or USA improved 40,000,000+ longer lifespan.
3. Safe, Fast, Painless, permanent
4. For all skin and hair style hair removal
5. 14.0 inch true color touch screen with 6 skin colors.
6. Two models can be selected, single handle or double handle



**Most advanced TEC cooling system inside this machine to deliver a super good cooling results. 24 hours working time.Sapphire material.Actual Pain-free experience!**

One minute cooling -21°C



## Best **TEC Cooling** system hair remover laser



- 1)10 bars/ 12 bars/ 16 bars, 10-40 million shots.
- 2)Big spot 12\*20mm / 12\*35mm for fast hair removal.
- 3)Frozen ice cooling sapphire crystal, for painless hair removal.
- 4)Hight quality gold laser bar imporved from German.
- 5)Medical handle as low as 350g easy to treatment.

**1 minute**  
and you're done

35mm {

4passive

**FASTER**

Screen system



1. 14 inch big touch screen, humanized designed.
2. Male/female, body parts, 6 skin types for option.
3. 6 language systems, native language + whatever you want.
4. Free OEM program service (logo/system...)

Laser bar	Germany/USA imported, 8 bar, 10 bar, 12 bars, 16 bars
Wavelength	755nm 808nm 1064nm/ or single 808nm
Energy	0-150J/cm <sup>2</sup>
Frequency	1-10hz
Pulse width	40-800ms
Operative LCD Interface	14.0 TFT True Color LCD touch screen
Cooling system	Air and Close Water Circulation Treatment head with semiconductor cooling temperature adjustable
Power	4000W
Electrical source	220V/50HZ or 110V/60HZ
N.W/ G.W	62kg/78kg

**Standardized  
production**



**Professional  
Production R & D  
Sales Department**



**Mass production  
of various styles**



- >>OEM your own machine body.
- >>OEM your own program software.
- >>OEM pretty machine colors.
- >>Print your LOGO on machine body and screen.

## Custom Section


**CASE** design



**LOGO** design



**SYSTEM** design



**Certificate Display**

# CERTIFICATE DISPLAY

