

The advertisement features a dark blue background with a futuristic, tunnel-like structure made of vertical, curved panels that recede into the distance. A large, glowing blue number '1' stands on a circular platform in the center. Above the '1', the text 'MLKJ LIFETIME SERVICE' is displayed in white, followed by 'FIRST COMPANY THAT OFFERED 2 YEARS WARRANTY' in a larger, bold white font. A blue button with the text 'START SERVICE' is positioned above the number '1'. In the top right corner, the MLKJ logo is visible. At the bottom, three service offerings are listed with icons: a wrench and screwdriver for '24-Hour online After-sales service', a factory icon for '15 Years Laser Beauty Device Manufacturer', and a smartphone with a checkmark for 'OEM Factories Support customization'.

MLKJ LIFETIME SERVICE

**FIRST COMPANY THAT OFFERED 2 YEARS WARRANTY**

START SERVICE

1

 **24-Hour online**  
After-sales service

 **15 Years Laser Beauty**  
Device Manufacturer

 **OEM Factories**  
Support customization

**First company that offered two years warranty□lifetime service.**

EXFU Seiko



# 11D HIFU

Macro & Micro High  
Intensity Focused Ultrasound System



SAFE



EFFECTIVE



MULTI-DEPTH



Skin rejuvenation laser body and face slimming machine skin rejuvenation system 11D Focused Ultrasound 11D HIFU

## Products Description

## ***Macro & Micro High Intensity Focused Ultrasound System***

Non-Surgical Treatment Of Face Lifting, Body Firming And  
Body Contouring System!



### SIZE OF CARTRIDGE

1.5mm,2.0mm,3.0mm,4.5mm for face treatment  
6.0mm.9mm.13mm for body treatment.



Use high energy ultrasound focusing on the treatment position, so that the skin tissue create the thermal, and make high-speed friction cells to stimulate collagen. Such heat effect will not hurt the epidermis, because the treatment fast and direct access to the treatment site within 0-0.5 seconds, without touching the surrounding tissue and can be passed directly to the superficial tendon membrane system (SMAS), so it can be done firm skin while pulling the muscle layer, thin face up to the progressive effect.

## **MF2.0MM Unique ultra-thin probe Clexis patented technology**


"Special cartridge for eyes treatment"



Accurate target focusing



# Body Treatment



## • *Treatment Effect* •



Name	High Intensity Focused Ultrasound
Handles	Facial: 1.5mm, 3.0mm, 4.5mm 11D HIFU: MF2/6/9/13 L4-4.5 L7-1.5 L7-3.0
Energy output	0.1J-2.5J (0.1J per step)
Space between Point	1.0mm-4.0mm
Shot length	5mm-25mm
shots	Each cartridge 300000 shots
Screen Size	15 inch
Voltage	110V-240V, 50/60Hz

## FEATURES



# SMART SYSTEM

GET TO WORK WITH ONE CLICK



Color touch screenskin rejuvenation

1. 15.6 inch big touch screen, humanized designed
2. Male/femal, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service(logo/system...)

**Welcome to consult more product details!**