



Exfu is committed to achieving efficient methods with low time costs.

Life is beautiful, some if there is a perfect solution, do it!

Maybe you don't have time to go to the gym, and the long process of getting your muscles into shape.

Then try our latest magnetic molding machine!

It has FDA professional approval and is even favored by fitness trainers!



EXFU FindX

In addition to the three best-in-class technologies, EXFU FindX is also equipped with additional features that can improve ease-of-use for operators and clinical outcomes for patients.

1. Four handles support continuous operation

2. One-click quick startup

3. Large 14-inch ultra-clear screen

Focused electromagnetic field with safe intensity to pass through the body and interact with motor neurons. Magnetic lean provides the body with severe

maximum contraction stimulation. The repeated contraction of muscle consumption metabolizes excess fat, reconstructs the internal structure, makes the production of muscle protein and muscle fiber chain, and creates higher density and muscle volume.

EXFU™ Find X Complete 360° Body Shaping Solution

Targeting Fat, Muscle,
and Body Treatments
The Advanced 5-in-1 Solution



5-in-1 independent work station, 4 high-strength applicators and 1 buttock augmentation cushion, fat reduction, shaping, lifting, buttock augmentation, all have good results.

Good result

30 minutes treatment= 30000 workouts, increase 16% muscle and reduce 21% fat after about 2-4 courses

7 Tesla High Intensity


cover big skeletal muscles, and remodel its inner structure.

4 +1 magnetic applicators

break down fat deposits and increase muscular tone and strength

Targeting 5 body parts

abdomen, arms, hips, hamstring, thighs



**ONE SESSION(30MINS/TREATMENT,
4 TREATMENT IN 2 WEEKS):**

1. 7 TESLA INTENSITY
2. 16% MUSCLE BUILDING/19% FAT REDUCTION
3. EQUIVALENT TO 20,000 SIT-UPS/50,000 MUSCLES CONTRACTIONS
4. DUAL SCULPTING FOR SAME/ DIFFERENT PERSON SAME TIME

Stronger Stimulation


Double layer Coil generates deeper and stronger high-intensity electromagnetic effects.

Liquid Cooling Technology

TeslaSculpt equipped with liquid-cooled applicators that allow long time operation without any overheating issue.

Four Sculpting

4 Magnetic Stimulation Applicators to treat all of your concerns independently or simultaneously.



EXFU FINDX

EXFU FindX modality offers an unprecedented solution to the muscle category. Controlled and precise energies are delivered uniquely to muscle groups via best-in-class electrical muscle stimulation. The muscles are safely activated without heating the surrounding tissue. Paired with specific pulse durations by muscle group, performance-informed smart training programs and advanced adaptive algorithms, each treatment is fully personalized and mimics a real training session that can build intensity over time. exfu enables effective and comfortable treatment programs regardless of the patients' fitness level, resulting in stronger and more defined muscle appearance.

5-in-1 independent work station, 4 high-strength applicators and 1 buttock augmentation cushion, fat reduction, shaping, lifting, buttock augmentation, all have good results.





With 4K HD screen, preset different treatment parameters for men and women, it can do personalized customized treatment data.

We also provide professional custom service, if you need specific language and operation technology, you can also contact us directly.



FDA Approved Ems Sculpting Machine 5 in 1 ems sculpt body slimming muscle stimulator machine