



C2 Vertical Magnetic Thin

Painless, sweat-free, convenient and fast

Multi-angle display



Simpler, higher smile
No diet, no exercise



How does this machine work?

It adopts high intensity pulsed electromagnetic to build muscle and tones for the abdomen, arms, legs, and buttocks.


Slimming your body and becoming more fit are easier than you think when you use the EMS beauty machine. This whole-body strengthening device can be used by men and women to quickly change their bodies. This machine allows you to tone your muscles and remove stubborn fat with a short session.

Bodybuilding fat burning


Get in shape easily




Simpler, higher smile
No diet, no exercise




Arm



Vest line



Peach hip




Beautiful legs



Muscle Building

The Intensive supramaximal contractions stimulate the react of the muscles. The energy reconstructs internal structure, causing the regrowth of muscle cells, the creation of protein, and muscle fiber chains. This procedure leads the results in higher density and muscle volume.



A variety of styles and specifications can be selected

Fat Removal

The muscle-building improves the epinephrine release, which brings an obvious effect that leads to supramaximal lipolysis in fat cells.

Circulation of lymphatic improvement

The electromagnetic waves stimulate the hemoglobin, increase blood flow, lymphatic circulation. The procedure increases that dead fat cells are eliminated naturally.

High Intensity


TeslaSculpt based on high-intensity magnetic energy, which could cover big skeletal muscles of the human body, and this high-energy level allows muscle response with a deep remodeling of its inner structure.

Working principle

Using non-invasive High-Intensity Focused Magnetic Vibration (HIFEM) technology, Magnetite uses two oversized treatment handles to release high-frequency magnetic vibration energy, which penetrates to a depth of 8cm muscle tissue, induces continuous expansion and contraction of muscles, and achieves high-frequency extreme training, deep remodeling and growth of myofibrils (muscle enlargement), and the production of new collagen chains and muscle fibers (muscle hyperplasia), thereby training and increasing muscle density and volume.

The 100% extreme muscle contraction of HIFEM technology can trigger a large amount of fat breakdown. The fatty acid concentration is too high, causing the fat cells to die and be excreted by the body's normal metabolism within a few weeks.

Therefore, Magnetic Lisu can strengthen and increase muscle, and at the same time achieve the effect of reducing fat.



Stronger Stimulation

Double layer Coil generates deeper and stronger high-intensity electromagnetic effects.

Liquid Cooling Technology

TeslaSculpt equipped with liquid-cooled applicators that allow long time operation without any overheating issue. 4 Magnetic Stimulation Applicators to treat all of your concerns independently or simultaneously.

Effect comparison



em slim machine
ems body sculpting machine
ems body slimming beauty machine
ems body slimming machine

**Please consult us for more product details
24 hours online!**