

NEW ARRIVAL



C2 Vertical Magnetic Thin

- ◆ BUILD MUSCLE
- ◆ BURN FAT

The newest air-cooling system
4 working handles





Advantages of weight loss machine

1. Lying Down for 30 Minutes=5.5 Hours of Exercise
2. Non-invasive treatment with no discomfort,Painless
3. 19% Reduction in Fat
4. 16% Increase in Muscle Mass
5. No downtime
6. Improve confidence in appearance
7. Noticeable, natural-looking results
8. Simple operation and Short treatment times



Operation



Hamstring



Thighs



Abdomen




Arm



Hip



BEFORE AND AFTER



Before

5 Weeks After

8 Weeks After



Before

After



Before

After







15 Years Laser Beauty Device Manufacturer

Make your beauty business easier



R & D design Design and development teams from the United States and Germany.

Certified and tested mature production line and quality inspection team and standard production line.

What we care about most is High-quality products are the standard for maintaining cooperation.

Our Factory Show

**15 Years
Professional Factory**




**Professional
Production R & D
Sales Department**



**Mass production
of various styles**



**Real Feedback
of Alibaba**



Certificate display

Certificate is the best guarantee for you



Packing And Shipping

Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service

Delivery time: 3-7 working days after payment

Delivery way: by air, DHL , TNT , Fedex, or delivery agent

Before delivery: testing video and pictures

After delivery Warranty: 2 years, offer impeccable after-sale service

CONTACT US



In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV SFDA ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month

FAQ

Q1. When will the machine be delivered to me?

A: We will make delivery within 7 days, after payment is made.

Q2. Do you provide training services?

A: We provide online training services, after training, we will send you a certificate.

Q3. What if the machine breaks down?

A: Under warranty, we will send you free components, the machine is module-designed.

Every part is easy to replace.

Q4. How to make payment?

A: We accept Alibaba Trade Assurance payment, wire transfer, Western Union, VISA.

Q5. What about guarantee?

A: 2 years for machine body, 1 year for handle. Under guarantee if the handle breaks down, we will send you a new handle.