




**Globally accessible tattoo removal machine portable picosecond tattoo removal machine with 5 heads tattoo removal system**

15 years beauty instrument production and sales manufacturers

Beauty devices have CE/FDA certification

24-48 hours delivery

Globally accessible



U.S. FOOD & DRUG

ADMINISTRATION

Follow FDA | En Español

SEARCH

Home

Food

Drugs

Medical Devices

Radiation-Emitting Products

Vaccines, Blood & Biologics

Animal & Veterinary

Cosmetics

Tobacco Products

Establishment Registration & Device Listing

FDA Home

Medical Devices

Databases

New Search

Back To Search Results

Establishment:

WEIFANG MINGLIANG ELECTRONICS CO.,LTD.

Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District

Weifang Shandong, CN 261000

Registration Number: 3007549439

FEI Number\*: 3007549439

Status: Active

Date Of Registration Status: 2022

Owner/Operator:

Weifang Mingliang Electronics CO.,LTD.

Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District

Weifang, Shandong CN 261000

Owner/Operator Number: 10082104

Official Correspondent:

Pingping Lv

Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District

Weifang, Shandong CN 261000

Phone: 86-10-56335780

US Agent:

Wen Xu

4611 Maydew Mews Way Ellicott City MD

Ellicott, MD US 21042

Phone: 443 8968583 Ext

Email: Infor@Believe-Med.Com

\* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English

Picosecond laser tattoo removal 532 755 1320 1064nm  
q switch nd yag laser tattoo removal machine

## Products Description

# Picosecond new upgrade

The maximum energy, the best effect

Smallest size


Rational design

5 treatment heads, detachable and rotatable outer wire design  
a solution to a variety of skin problems



## Craftsman quality, Reflect the extreme details Rational design

1. The product has a novel appearance and is mainly sold as a hot model
2. Adjustable handle bracket, strong and durable
3. Non-slip handle, ergonomic design
4. Extended treatment connecting tube



### **Smallest size**

1. Adopt core components of medical equipment. The size of the machine is smaller.
2. Lower cross-border freight than ordinary machines.
3. The handle is smaller and more convenient to use

### **About Handles**

## How is it different?



### **5 treatment heads, detachable and rotatable outer wire design**

1. 532: Remove colored tattoos.
2. 755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin
3. Adjustable parameter 1064: for tattoo removal, eyebrow clearing
4. Non-adjustable parameter 1064: used for removing tattoos and clearing eyebrows
5. 1320: large area treatment, large spot, fast treatment speed and good effect



## Advantage

### The maximum energy, the best effect


Laboratory data shows that the laser rod reaches 2000mj.  
Precise calculation of laser value  
The aiming distance of the indicator light is 6-10cm

### 5 treatment heads

Detachable and rotatable outer wire design  
532: Remove colored tattoos.  
755:honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin  
Adjustable parameter 1064: for tattoo removal, eyebrow clearing  
1064: used for removing tattoos and clearing eyebrows  
1320: large area treatment, large spot, fast treatment speed and good effect

## Effect comparison

# Working Principle



## TREATMENT HEADS



532: Remove colored tattoos.



755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin



Adjustable parameter 1064: for tattoo removal, eyebrow clearing



1064: used for removing tattoos and clearing eyebrows



1320: large area treatment, large spot, fast treatment speed and good effect

## Product Features



Tattoo removal

Eyebrow embroidery removal

Freline removal

Eyeline removal

Whitening

Skin rejuvenation

Lipline removal

Pigment treatment

### The maximum energy, the best effect

1. Laboratory data shows that the laser rod reaches 2000mj.
2. Precise calculation of laser value
3. The aiming distance of the indicator light is 6-10cm


### Product Paramenters

1064nm

## Precise Calculation Laser Value

The aiming distance of the indicator light is 6-10cm

6-10cm



You think you can. you can  
No one can stop me





Imported stone Q-switch cassette mechanism, the entire solid laser, without replacing Q-switch.

Tattoo removal. Remove black, blue, brown and red tattoo pigments on eyebrow, eye line, lip line and other parts of body Remove birth mark, otas nevus and dermal speckle  
Instantaneous blasting high laser technology, international standard of production, technical test strictly.

A solutionto avariety of skin problems.

# Treatment **Effec**

Tattoo removal

Nevus Ota removal

Eyebrow and eyeliner removal

Lipline removal



## Screen Display

10 inch true color big touch screen


Super cool touch Experience

6 language for your choice, Special language can be added

OEM service are accepted,such as adjust parameter,include energy/fre/time scope

# Screen Image

The display screen is clear at a glance and easy to operate



**10** inch true color big touch screen. Super cool touch Experience.  
Special language can be added according to your actual requirements.  
Your logo also can be added on diode laser screen.



# Picosecond new upgrade

The maximum energy, the best effect

Smallest size

Rational design

5 treatment heads, detachable and rotatable outer wire design  
a solution to a variety of skin problems



Technology	Pico second laser / Nd yag laser
Function	Tattoo freckle acne removal / Blackface doll treatment
Handles	5 Heads
Wavelength	1064nm 532nm 1320nm 755nm
Treatment area	Face Body Face Neck Arms Legs....
Frequency	1-10HZ
Voltage	110V/220V 50-60Hz

## Why choose us

15 years manufacturer with rich experience in producing beauty machines!

15 years Alibaba golden supplier

CE certificates approved TUV certificates approved Personal custom made & OEM





**15 Years**  
Laser Beauty Device  
**Manufacturer**

Experienced manufacturer You will cooperate with a factory with more than 15 years of experience in manufacturing beauty equipment, integrating R&D, design, production and sales. We exported to 126+ countries, supporting 60,000+ users, more reliable and more professional.

Professional R&D and newest technology Relying on our professional R&D and technical team, products will continue to be upgraded, you will get the newest machines with latest technologies, and always be at the forefront of the industry, making you more competitive in the beauty industry.

# Our Showroom



All kinds of certifications All the beauty machines you get have passed high quality certifications, such as CE, FDA, TUV, ISO 13485...etc.

No matter which country/area you are from, you can use the machines legally, safely and effectively.

Excellent ODM/OEM service You will receive free customization services and exclusive protection to provide a strong guarantee for your brand development. Your machines can be unique on the market.

# CERTIFICATES



BV-France's most authoritative  
certification body BV Group certification



TUV-Germany's most authoritative  
certification body Rheinland certification



SGS-Switzerland's most authoritative  
certification body SGS certification



Italy ECM Certification



# ISO 13485

NQA medical system  
certification





# Production Process



Warehouse spare parts



Parts assembly



Precision machining



Handle fabrication



Product assembly process



Quality inspection  
of finished products



Product packaging

No worry about warranty and after sale service All machines have a two-year warranty and lifetime maintenance, you get 24-hour online service from a professional after-sales team. Whenever you get any problems, you can contact our after-sales team, and the problem will be solved in the fastest time.

Free training on all machines You will receive free training and training certificate. Get complete and easy-to-understand training courses. Also professional beauticians will provide online guidance. Your beauty business will become easier.

## Product packaging



# Packing

Safe and thoughtful packaging service



## Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board  
inside aluminum box**



## 3-5 business days delivery

Professional guide:

24 Hours Professional online customer service

# CONTACT US



## In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV SFDA ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month