

Already help over customers from 102 countries make good business.

Various templates for you to choose.

Product Description




Advantages


- 1.3 waves Yag 1064nm +Diode 808nm + Alex 755nm, 3 in 1 all skin type can be used
High power: NEW 1000W/ 1200W /1600W handle for options
- 2.Big spot size 15*25mm, 12*20mm, 12*38mm, face & body hair removal solution!
- 3.Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
- 4.Long life: double life support guarantee



A natural



B energy absorbed by hair



C hair follicle died

Safe Use On All Skin Type



Alex 755nm for light skin

Speed Diode 808nm for yellow skin

Long Pulse 1064nm for darker skin

Handle Details

Three Wavelengths Work Together




Features & Diode Laser

1. Laser bar: Imported laser bar US Coherent/ Dilas brand
2. 3 waves: 755nm 808nm 1064nm 3 in 1
3. Spot size: 12*20mm, 12*35mm
4. Power: High power 1000W 1200W 1600W & OEM
3. Cooling: Super cooling start frozen in 1 minutes, 24 hr continue work

Software & Program

1. Easy Start: Based on male & femal/6 body parts, Add 6 skin types setting smart software
2. 6 language options, your native language can be added
3. OEM program service available(logo/system...)





**For more model details, please contact us!
Online 24 hours!!**