



Product Advantage: skin rejuvenation

*IPL handle use UK lamp, 150,000-30,000 shots. Spot size 15*50mm. Cooling performance is excellent.

*It use 15.6 HD screen super sensitive color LCD touch screen, all parameters setup can be finished easily here.

*Add your own logo and language into the machine.

Products Description

IPL+Laser

Multifunctional Machine

IPL handle use UK lamp, 150,000-30,000 shots.
Spot size 15*50mm. Cooling performance is excellent.



| 420nm | 510nm | 560nm | 585nm | 640nm | 690nm | 755nm |
|------------------------------------|---------------------------------------|--------------------------------------|-----------------------|---------------------------|---------------------------|---------------------------|
| Acne removal Vascular treatment | Pigmentation therapy, Remove freckles | Skin rejuvenation wrinkle removal | Remove Red blood silk | Hair removal (light skin) | Hair removal (brown skin) | Hair removal (black skin) |
| | | | | | | |

Color touch screenskin rejuvenation,One-button startup

1. 15.6 inch big touch screen, humanized designed
2. Male/femal, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service(logo/system...)

SMART SYSTEM

Full screen, Wide field of view

15.6" HD SCREEN



EXFU

One touch power on

ipl hair removal and skin rejuvenation device
ipl hair removal and skin rejuvenation new
ipl hair removal skin rejuvenation machine



ipl hair removal and skin rejuvenation machine

ipl hair removal skin rejuvenation machine

ipl rejuvenation machine

24-hour online service!Welcome to consult!