



Product Advantage: skin rejuvenation

IPL handle use UK lamp, 150,000-30,000 shots. Spot size 15*50mm. Cooling performance is excellent. It use 8.4TFT super sensitive color LCD touch screen, all parameters setup can be finished easily here. Add your own logo and language into the machine.



About Handle

Elight Handle Design



Specification

Optional Spot Size:
15*50mm

Imported UK lamp

Long shots time:
150,000--300,000 times

7filters:
420/510/560/585/640/690/755nm

Super cooling system



Suitable to all types of skins ,high efficient treatmentresults.



8.4 inch color touch screenskin rejuvenation

1. 8.4inch big touch screen, humanized designed

2. Male/female, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service(logo/system...)

Screen interfac

Professional Screen, Easy operation



Main Menu

8.4 inch true color big touch screen
Super cool touch Experience.



IPL mode:

5 options for skin treatment



IPL mode:

operation interface




Multi-language screen:

6 languages for your choice



System setting



System setting

Size Information

Probe	IPL: 15*50mm imported UK lamp
Wavelength	420/510/560/585/640/690/755nm
Lifetime	150,000-30,000 shots
Screen	8.4 inch Colorful Touch Screen
Cooling System	1200W TEC Cooling + Air Cooling + Water Cooling 3 in 1 system
Warranty	2 Years for machine
Plug	US,EU,AU, 110V and 220V
Weight	30KG
Logo	Add logo for free
Language	Any language as yours

Multi-angle display



Left view



Front view



right elevation