

Various Modes Of Operation

Machine, mobile phone, pad,
synchronization control



Mobile Phone



Pad



Machine



12*24mm or 12*36mm handle

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...



Advantages

- 1.3 waves 1064nm + 808nm + 755nm, 3 in 1 all skin type can be used high power: NEW 1200W /1800W handle for options
 - 2.Big spot size 12*20mm, 12*36mm, face & body hair removal solution!
 - 3.Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
 - 4.Long life: 40 million shots
- Customized service for all style all colors as you want


Function of diode laser

The 3 waves works in the same handle at the same time

755nm for white skin(fine, golden hair)

808nm for yellow/neutral skin

1064nm for black(black hair)



Super ice cooling effect, painless treatment is more comfortable



Machine Display




Screen Display

1. 14 inch big touch screen, humanized designed
2. Male/female, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service (logo/system...)

The new **14-inch 4K** screen


Support system customization



Detailed specifications of diode laser ice platinum

Item	Vertical 755nm+808nm+1064nm diode laser
Laser type	Diode laser,10bars,12bars,16bars
Power of each bar	120w
Display system	14 inch true color touch screen
Wavelength	808nm/755nm/1064nm
Spot Size	12*24mm ² 12*36mm ²
Laser handle lifespan	Max 40,000,000 shots
Energy	10-150 J / cm ² adjustable (up to 150J/cm ² optional)
Pulse Duration	10~400ms adjustable
Frequency	1.0-10 Hz;
Cooling system	TEC semiconductor condenser+air+close water circulation
Handle Power	1000w/1200w/1600w
Machine Power	2000w
Electrical source	220V/50HZ or 110V/60HZ
Crystal temperature	-20~5°C
Machine temperature	24°C constant temperature
Power requirement	AC220V/110V 50~60Hz
Gross weight	80KG

Our factory


**U.S. FOOD & DRUG
ADMINISTRATION**


Follow FDA | En Español

SEARCH

[Home](#) | [Food](#) | [Drugs](#) | [Medical Devices](#) | [Radiation-Emitting Products](#) | [Vaccines, Blood & Biologics](#) | [Animal & Veterinary](#) | [Cosmetics](#) | [Tobacco Products](#)

510(k) Premarket Notification

[FDA Home](#) | [Medical Devices](#) | [Databases](#)



[510\(k\)](#) | [DeNovo](#) | [Registration & Listing](#) | [Adverse Events](#) | [Recalls](#) | [PMA](#) | [HDE](#) | [Classification](#) | [Standards](#)
[CFR Title 21](#) | [Radiation-Emitting Products](#) | [X-Ray Assembler](#) | [Medsun Reports](#) | [CLIA](#) | [TPLC](#)

[New Search](#)
[Back To Search Results](#)

Device Classification Name	Powered Laser Surgical Instrument
510(K) Number	K200525
Device Name	Medical Diode Laser Hair Removal System
Applicant	Weifang Mingliang Electronics CO., LTD No. 13426, Yuqing Street, High-Tech District Weifang, CN 261000
Applicant Contact	Li Mingliang
Correspondent	Beijing Believe-Med Technology Service Co., Ltd. Rm. 912, Building#15, XiYueHui, No.5, YiHe North Rd., FangShan District Beijing, CN 102401
Correspondent Contact	Ray Wang
Regulation Number	878.4810
Classification Product Code	GEX
Date Received	03/02/2020
Decision Date	05/27/2020
Decision	Substantially Equivalent (SESE)
Regulation Medical Specialty	General & Plastic Surgery
510k Review Panel	General & Plastic Surgery
Type	Traditional
Reviewed By Third Party	No
Combination Product	No

Custom Section

CASE design



LOGO design



SYSTEM design

