

# Soprano ICE PLATINUM

## Medical beauty machine

Medical diode laser hair removal  
For all skin types,  
triple wavelength laser work together

### Diode laser hair removal machine.

- \* No pigmentation. Suitable for any kind of skin;
- \* Suitable for any kind of skin ;
- \* Safe and Fast.
- \* International Hair Removal Golden Standard;
- \* The optimal laser wavelength 808nm or 755+808+1064nm;

**350g**

Ultra light handle



12×36mm

**Handle**

1800W diode laser hair removal




12×24 mm

**Handle**

1200W diode laser hair removal

**Three Waves laser combines 755nm 808nm 1064nm, suits for all skin and hair types treatment.**

755nm for white skin (fine/ golden hair) hair removal.  
 808nm for yellow (natural) skin hair removal.  
 1064nm for black skin (black hair) hair removal.  
 Before and After ice platinum diode laser hair removal machine




**ALEX 755nm**  
**For white skin**  
 (fine, golden hair removal )



**SPEED 810nm**  
**For yellow neutral skin**  
 (brown hair removal)



**YAG 1064nm**  
**For black skin**  
 (black hair removal)



# Treatment Comparison

