



Medical

1200W+1800W high power handle
14-inch 4K HD large screen
755nm+808nm+1064nm
Painless permanent hair removal

Diode Laser

New medical diode laser hair removal machine

14-inch 4K HD large screen

1200W+1800W high power medical handle

Handle design is more technological

The handle weighs about 350g, which is smaller and lighter

12 * 24mm + 12 * 36mm spot, fully meet the treatment needs of different parts

Machine Display



- 1.Stronger, faster, smarter and most importantly pain free approach for hair removal.
2. 755nm, 808nm, 1064nm Diode Laser hair removal machine
- 3.TEC cooling system,1 minute handle cooling



Before and After

Treatment comparison



100% effective device without overriding skin tone, sugar, hair color. For more application details, please consult !