



Medical

1200W+1800W high power handle  
14-inch 4K HD large screen  
755nm+808nm+1064nm  
Painless permanent hair removal

Diode Laser

New medical diode laser hair removal machine

14-inch 4K HD large screen

1200W+1800W high power medical handle

Handles design is more technological

The handle weighs about 350g, which is smaller and lighter

12 \* 24mm + 12 \* 36mm spot, fully meet the treatment needs of different parts

# Machine Display



- 1.Stronger, faster, smarter and most importantly pain free approach for hair removal.
2. 755nm, 808nm, 1064nm Diode Laser hair removal machine
- 3.TEC cooling system,1 minute handle cooling



# Before and After

Treatment comparison



100% effective device without overriding skin tone, sugar, hair color. For more application details, please consult !