



Hair Removal + ND YAG

*2in1 Multifunctional combination machine
one machine and two effects*

Multi-functional beauty instrument manufacturers

2 in 1 Diode laser hair removal +nd yag laser tattoo removal

Products Description

2in1 Multifunctional combination machine




- ✓ Hair removal
- ✓ Toner whitening
- ✓ Mole removal



- ✓ Tattoo removal
- ✓ Freckle
- ✓ Birthmark removal

Real shot display

Higher-definition screen, stronger cooling effect,
more effective treatment



Diode laser+Elight/IPL+Nd yag+RF, 4 different handles into one machine

1.This 4 in 1 RF+Elight IPL+diode laser+nd yag laser have high quality accessories , strong power*Diode laser handle use Germanyimported golden laser bar ,50 million shots time.


*Elight handle use UK lamp, 150,000-30,000shots

*Nd yag :Equip 5 treatment heads, adjustable 1064nm/532nm, fixed 1064nm/532nm, fixed 1320nm.


*RF :Equip 6 treatment heads, used for wrinkle removal, skin lifting, skin tightening.

Three wavelengths work together

12*24mm² {



12*35mm² {



Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.





HANDLE-ADVANTAGE



specification

- **Multiple choice**
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- **Multifunction**
1320nm for skin rejuvenation
- **Indicator light**
Dot position is the indicator light position dot more accurate



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm



Fixed
1064nm



Before and after



14inch true touch screen

6 kinds skin type for choice. accept ODM OEM Service.

The new **14-inch** screen is clear easy to operate.



For more information about product introduction and use details,

please consult us immediately

Online 24 hours!!