



WORK THEORY

1 Non-invasive HIEMT procedure induces powerful muscle contractions not achievable through voluntary contractions. When exposed to strong contractions, the muscle tissue is forced to adapt to such extreme condition. It responds with a deep remodeling of its inner structure that results in muscle building and sculpting your body.

2 At the same time, the 100% extreme muscle contraction of HI-EMT technology can trigger a large amount of fat decomposition excreted by the body's normal metabolism within a few weeks. Therefore, slim beauty machine can strengthen and increase muscle, and reduce fat at the same time.

The image shows the EXFU FindX machine, a sleek, modern device with a large 14-inch screen and two handles. The machine is white with blue accents. To the right of the machine, the text "EXFU FindX" is displayed in a large, bold font. Below the title, there is a paragraph of text and a list of three features. The background is dark with a pattern of curved lines.

EXFU FindX

In addition to the three best-in-class technologies, EXFU FindX is also equipped with additional features that can improve ease-of-use for operators and clinical outcomes for patients.

1. Four handles support continuous operation
2. One-click quick startup
3. Large 14-inch ultra-clear screen

Focused electromagnetic field with safe intensity to pass through the body and interact with motor neurons. Magnetic lean provides the body with severe maximum contraction stimulation. The repeated contraction of muscle consumption metabolizes excess fat, reconstructs the internal structure, makes the production of muscle protein and muscle fiber chain, and creates higher density and muscle volume.

EXFU™ Find X Complete 360° Body Shaping Solution

Targeting Fat, Muscle,
and Body Treatments
The Advanced 5-in-1 Solution



5-in-1 independent work station, 4 high-strength applicators and 1 buttock augmentation cushion, fat reduction, shaping, lifting, buttock augmentation, all have good results.

Good result

30 minutes treatment= 30000 workouts, increase 16% muscle and reduce 21% fat after about 2-4 courses

7 Tesla High Intensity

cover big skeletal muscles, and remodel its inner structure.

4 +1 magnetic applicators

break down fat deposits and increase muscular tone and strength

Targeting 5 body parts

abdomen, arms, hips, hamstring, thighs



Stronger Stimulation


Double layer Coil generates deeper and stronger high-intensity electromagnetic effects.

Liquid Cooling Technology

TeslaSculpt equipped with liquid-cooled applicators that allow long time operation without any overheating issue.

Four Sculpting

4 Magnetic Stimulation Applicators to treat all of your concerns independently or simultaneously.



5-in-1 independent work station, 4 high-strength applicators and 1 buttock augmentation cushion, fat reduction, shaping, lifting, buttock augmentation, all have good results.





With 4K HD screen, preset different treatment parameters for men and women, it can do personalized customized treatment data.

We also provide professional custom service, if you need specific language and operation technology, you can also contact us directly.



FDA Approved Ems Sculpting Machine 5 in 1 ems sculpt body slimming muscle stimulator machine