




Ems slimming machine hiemt pro machine Build muscle build muscle ems body sculpting machineem for sale



How does this machine work?

It adopts high intensity pulsed electromagnetic to build muscle and tones for the abdomen, arms, legs, and buttocks




Four handles work together

30 Mins Treatment
= 30,000 Workouts

Muscle Building

The Intensive supramaximal contractions stimulate the react of the muscles. The energy reconstructs internal structure, causing the regrowth of muscle cells, the creation of protein, and muscle fiber chains. This procedure leads the results in higher density and muscle volume.



A variety of styles and specifications can be selected

Fat Removal

The muscle-building improves the epinephrine release, which brings an obvious effect that leads to supramaximal lipolysis in fat cells.

Circulation of lymphatic improvement

The electromagnetic waves stimulate the hemoglobin, increase blood flow, lymphatic circulation. The procedure increases that dead fat cells are eliminated naturally.

High Intensity

TeslaSculpt based on high-intensity magnetic energy, which could cover big skeletal muscles of the human body, and this high-energy level allows muscle response with a deep remodeling of its inner structure.

SCREEN
4K HD



Stronger Stimulation

Double layer Coil generates deeper and stronger high-intensity electromagnetic effects

Liquid Cooling Technology

TeslaSculpt equipped with liquid-cooled applicators that allow long time operation without any overheating issue

4 Magnetic Stimulation Applicators to treat all of your concerns independently or simultaneously

Effect comparison




Before

Five weeks later

Eight weeks later



Before



After a month



Before



After a month



Before



After a month



Before



After a month

**Please consult us for more product details
24 hours online!**