

# Various Modes Of Operation

Machine, mobile phone, pad,  
synchronization control



Mobile Phone



Pad



Machine



# 12\*24mm or 12\*36mm handle

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...



## Advantages

- 1. 3 waves 1064nm + 808nm + 755nm, 3 in 1 all skin type can be used high power: NEW 1200W /1800W handle for options
  - 2. Big spot size 12\*20mm, 12\*36mm, face & body hair removal solution!
  - 3. Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
  - 4. Long life: 40 million shots
- Customized service for all style all colors as you want

## 755&808&1064nm Diode laser handle:

- \* Germany Dilas imported laser bar. (Newest and Best Quality Golden Laser Bar)
  - \* Shots: 10-40 million shots time
  - \* 755 808 1064nm is for all skin and hair style hair removal
- 1200w, 1600w, 1800W handle  
Same technology with the ice platinum




### Function of diode laser:

The 3 waves work in the same handle at the same time

755nm for white skin( fine, golden hair)


808nm for yellow/neutral skin

1064nm for black(black hair)



**Super ice cooling effect, painless treatment is more comfortable**





## Machine Display




### Screen Display

1. 14 inch big touch screen, humanized designed
2. Male/female, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service (logo/system...)

# The new **14-inch 4K** screen


## Support system customization



### Detailed specifications of diode laser

Item	Vertical 755nm+808nm+1064nm diode laser
Laser type	Diode laser,10bars,12bars,16bars
Power of each bar	120w
Display system	14 inch true color touch screen
Wavelength	808nm/755nm/1064nm
Spot Size	12*24mm <sup>2</sup> 12*36mm <sup>2</sup>
Laser handle lifespan	Max 40,000,000 shots
Energy	10-150 J / cm <sup>2</sup> adjustable (up to 150J/cm <sup>2</sup> optional)
Pulse Duration	10~400ms adjustable
Frequency	1.0-10 Hz;
Cooling system	TEC semiconductor condenser+air+close water circulation
Handle Power	1000w/1200w/1600w
Machine Power	2000w
Electrical source	220V/50HZ or 110V/60HZ
Crystal temperature	-20~5°C
Machine temperature	24°C constant temperature
Power requirement	AC220V/110V 50~60Hz
Gross weight	80KG

## Our factory

**U.S. FOOD & DRUG  
ADMINISTRATION**


Follow FDA | En Español

SEARCH

Home | Food | Drugs | Medical Devices | Radiation-Emitting Products | Vaccines, Blood & Biologics | Animal & Veterinary | Cosmetics | Tobacco Products

### 510(k) Premarket Notification

FDA Home | Medical Devices | Databases

510(k) | DeNovo | Registration & Listing | Adverse Events | Recalls | PMA | HDE | Classification | Standards  
CFR Title 21 | Radiation-Emitting Products | X-Ray Assembler | Medsun Reports | CLIA | TPLC

[New Search](#)[Back To Search Results](#)

Device Classification Name	<a href="#">Powered Laser Surgical Instrument</a>
510(K) Number	K200525
Device Name	Medical Diode Laser Hair Removal System
<b>Applicant</b>	<b>Weifang Mingliang Electronics CO., LTD</b>
	No.13426, Yuqing Street, High-Tech District Weifang, CN 261000
Applicant Contact	Li Mingliang
Correspondent	Beijing Believe-Med Technology Service Co., Ltd. Rm. 912, Building#15, XiYueHui, No.5, YiHe North Rd., FangShan District Beijing, CN 102401
Correspondent Contact	Ray Wang
Regulation Number	<a href="#">878.4810</a>
Classification Product Code	<a href="#">GEX</a>
Date Received	03/02/2020
Decision Date	05/27/2020
Decision	Substantially Equivalent (SESE)
Regulation Medical Specialty	General & Plastic Surgery
510k Review Panel	General & Plastic Surgery
Type	Traditional
Reviewed By Third Party	No
Combination Product	No

## Custom Section

### CASE design



### LOGO design



### SYSTEM design

