

Product Description

Advantages



1.3 waves 1064nm +Diode 808nm + 755nm, 3 in 1 all skin type can be used

High power: NEW 1000W/ 1200W /1600W handle for options

2.Big spot size 12*20mm, 12*35mm face & body hair removal solution!

3.Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.

4.Long life: double life support guarantee

Handle of diode laser hair removal machine:



Different spot size Handle color optional



Safe Use On All Skin Type


Alex 755nm for light skin
Speed Diode 808nm for yellow skin
Long Pulse 1064nm for darker skin

Treatment effect



Software & Program

1. Male & female, 6 body parts, 6 skin types setting smart software
2. 6 language options, your native language can be added
3. OEM program service available (logo/system...)



Product name	Laser Ice platinum hair removal machine
Laser bar	Germany/USA imported, 8 bar, 10 bar, 12 bars, 16 bars
Wavelength	755nm 808nm 1064nm/ or single 808nm
Energy	0-150J/cm2
Frequency	1-10hz
Pulse width	