

First company that offered two years warranty & lifetime service.

Body And Face Slimming Machine Wrinkle Remover Skin Revitalizer 11D Focused Ultrasound Machine Body And Face Slimming Machine 11d Focused Ultrasound





Skin rejuvenation laser body and face slimming machine skin rejuvenation system 11D Focused Ultrasound 11D HIFU

*A high repetition rate of short pulses are delivered deep into the dermis, achieving high average power and therapeutically effective heat build-up, with virtually no pain. *The latest machine has total 11 probes, which are fully carved and anti-aging. *Non surgical treatment of facial lifting, tightening and

body contour shaping system. The world's latest high energy focused ultrasound technology.

Products Description






MF2.0MM

Unique ultra-thin
probe Clexis patented
technology

"Special cartridge for eyes treatment"



Use high energy ultrasound focusing on the treatment position, so that the skin tissue create the thermal, and make high-speed friction cells to stimulate collagen. Such heat effect will not hurt the epidermis, because the treatment fast and direct access to the treatment site within 0-0.5 seconds, without touching the surrounding tissue and can be passed directly to the superficial tendon membrane system (SMAS), so it can be done firm skin while pulling the muscle layer, thin face up to the progressive effect.




Macro & Micro High Intensity Focused Ultrasound System

7CARTRIDGES+4 GUN HEADS

1.5mm,2.0mm,3.0mm,4.5mm for face treatment
6.0mm.9mm.13mm for body treatment.

Body Treatment



• *Treatment Effect* •



Model	7D Therapy Device Focused Ultrasound
Name	High Intensity Focused Ultrasound

Handles	Facial: 1.5mm, 3.0mm, 4.5mm
	7D HF: MF2/6/9/13
	L4-4.5 L7-1.5 L7-3.0
Energy output	0.1J-2.5J (0.1J per step)
Space between Point	1.0mm-4.0mm
Shot length	5mm-25mm
shots	Each cartridge 300000 shots
Screen Size	15 inch
Voltage	110V-240V, 50/60Hz



Color touch screenskin rejuvenation

1. 15.6 inch big touch screen, humanized designed
2. Male/femal, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service(logo/system...)

High-Intensity Focused Ultrasound

Anti-Wrinkle & V- shape face rebuilt



If you have any questions about beauty equipment, looking for partners and beauty entrepreneurship, please consult MLKJ!

Tel/Whatsapp: +86 15965364558

E-mail: sale@mingliangkj.com