



## Advantage

- ROHS & Medical CE Certified Machine
- 3 in 1 - Three Laser Technologies in one machine
- Alexandrite, diode and Nd:YAG lasers - treats all skin types
- Easy to operate
- Online training available
- 50 Million shots per handpiece!

Handle display

Ergonomic design

12°24mm

755nm 808nm 1064nm

Comfortable  
painless hair removal

Mild cooling  
No irritation

+

Optional flare

+

Cooling -25°/-23°C

ice

ONLY 350G  
HANDLE

# Medical diode laser

**50000000+shots**  
**can treat 10000+clients**

- \*Not easy to be burned
- \*Low maintenance cost
- \*Stable quality



\* Monitors applicator temperature 1,000 times per second

\*Maintains temperature throughout the session to deliver a comfortable, safe treatment. \*No downtime for cooling translates to higher productivity for your clinic

\*POWERFUL & EFFICIENT SYSTEM FOR ALL SKIN TYPES \*Optimal wavelength 755nm 808nm 1064nm

\*High peak power of 2000W,Fluence of up to 150 J/cm<sup>2</sup>/cm<sup>2</sup> \*Spot size 12x24mm

Screen system



\*Large touch color screen with advanced intuitive and easy human-computer interaction software.

\*Add your language for free.

Treatment effect





# *Before and After*

## Treatment comparison



- 1.Pain-free.The tip of handpiece cools down to 0-3C,makes treatment comfortable and effective.
- 2.Early treatment results.2-3 treatment sessions can show permanent results
- 3.Virtually pain-free contacted with sapphire of Real-5C more incomparable comfortable,and safety.

Company Profile





High Feedback Top10 Supplier  
OEM & ODM Professional Design



### Authentication

NQA13485,CB certification  
and more Complete certification  
materials



### 60000+ Customers

We have professional after sale  
team and we provide 7\*24  
services



### 2 Years Warranty

2 years warranty for each  
machine, All questions will get  
reponse within 24 hours

Weifang Mingliang Electronics Co., Ltd (abbreviation:MLKJ) located in Shandong, China. We have been in this field for more than 15 years, have tons of experience in design and manufacturing. Our most popular machines including 808nm diode laser, IPL machine, SHR machine, Nd yag laser, cavitation body slimming machine, picosecond laser... have been sold to more than 102 countries. We're so delighted that our products has brought value to our clients in the world.

***Please send me an inquiry for the BEST price!***