



Products Description





Product Advantage: skin rejuvenation

*IPL handle use UK lamp, 150,000-30,000 shots. Spot size 15*50mm. Cooling performance is excellent. It use 15.6 HD screen super sensitive color LCD touch screen, all parameters setup can be finished easily here.

*Add your own logo and language into the machine.



Multi-angle display



Non-cooling effect



Cooling effect

Specification

- High-quality IPL handpiece with super cooling system
- Optional Spot Size 15*50mm, 15*35mm, 12*30mm
- Long shots time:
IPL: 150,000 times
- 7 filters: 420/510/560/585/640/690/755nm
Standard: 420/510/560/690nm

Cooling Effect

The handle adopts a refrigeration effect which can effectively cool the skin in the treatment process and improve the comfort of treatments.


large
area

Cool
painless
comfortable








Protect
surrounding
skin

Fast
hair
removal

Filters for choice



ABS plastic filter box,
more exquisite, more durable

						
420nm	510nm	560nm	585nm	640nm	690nm	755nm
Acne removal Vascular treatment	Freckles removal	Skin rejuvenation wrinkle removal	Removal of visible blood vessels	Hair removal (fair skin)	Hair removal (medium skin tone)	Hair removal (darker skin)

Standard configuration: 420nm 510nm 560nm 690nm

Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

SMART SYSTEM

Full screen, Wide field of view

15.6" HD SCREEN



Color touch screenskin rejuvenation

1. 15.6 inch big touch screen, humanized designed
2. Male/female, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service (logo/system...)

24-hour online service!

Welcome to consult!