

## Exceeding your expectations

The most professional beauty machine customization service

## Surrounded by a whole set of services

Making the purchasing experience wonderful



**Surrounded by a whole set of services, making the purchasing experience wonderful**

**Exceeding your expectations, the most professional beauty machine customization service**


## High-Intensity Focused Ultrasound

Anti-Wrinkle & V- shape face rebuilt



**Skin rejuvenation laser body and face slimming machine skin rejuvenation system 11D Focused Ultrasound 11D HIFU**

## Products Description



**Macro & Micro High Intensity Focused Ultrasound System**

**7CARTRIDGES+4 GUN HEADS**


1.5mm,2.0mm,3.0mm,4.5mm for face treatment  
6.0mm.9mm.13mm for body treatment.




## MF2.0MM

Unique ultra-thin  
probe Clexis patented  
technology


"Special cartridge for eyes treatment"



Use high energy ultrasound focusing on the treatment position, so that the skin tissue create the thermal, and make high-speed friction cells to stimulate collagen. Such heat effect will not hurt the epidermis, because the treatment fast and direct access to the treatment site within 0-0.5 seconds, without touching the surrounding tissue and can be passed directly to the superficial tendon membrane system (SMAS), so it can be done firm skin while pulling the muscle layer, thin face up to the progressive effect.



# Body Treatment



## • *Treatment Effect* •



Name	High Intensity Focused Ultrasound
Handles	Facial: 1.5mm, 3.0mm, 4.5mm
	7D HF: MF2/6/9/13 L4-4.5 L7-1.5 L7-3.0
Energy output	0.1J-2.5J (0.1J per step)
Space between Point	1.0mm-4.0mm
Shot length	5mm-25mm
shots	Each cartridge 300000 shots
Screen Size	15 inch
Voltage	110V-240V, 50/60Hz

## ***Comparison before and after use***



Color touch screenskin rejuvenation

1. 15.6 inch big touch screen, humanized designed
2. Male/femal, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service(logo/system...)

EXFU Seiko



# 11D HIFU

Macro & Micro High Intensity Focused Ultrasound System



SAFE



EFFECTIVE



MULTI-DEPTH



Welcome to consult more product details!