



Diode laser+Nd yag laser

- 1.This 2 in 1 diode laser+nd yag laser have high quality accessories , strong power
- 2.Shots time for handles : Diode laser handle use Germany imported golden laser bar ,10-40 million shots time.
- 3.Nd yag :Equip 5 treatment heads, adjustable 1064nm/532nm, fixed 1064nm/532nm, fixed 1320nm.







About Handle

Three wavelengths work together

12*20mm² {



12*35mm² {



**Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.**



755+808+1064 diode laser suits for all skin and hair types treatment.

Imported Laser Bar. 20-40 million shots time.



Before

After




Before



After



Before



After

5 therapeutic **Heads**

■ 532 & 1064: removing tattoo, eyeliner and lip line

■ 1320: black-faced doll whitening and whitening skin

■ Adjustable 532 & 1064: removing large-area tattoos



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm




Fixed
1064nm

Nd yag tattoo removal :

Equip 5 treatment heads, adjustable 1064nm/532nm, fixed 1064nm/532nm, fixed 1320nm.

Working Principle






1. Male/femal, body parts, 6 skin types for option.

2. English+your native language.

3. Free OEM program service(logo/system...)

System Design



808 Mode:Operation interface



Laser Mode:Tattoo removal

Size Information

Item	808+nd yag laser
Handpiece	Nd yag laser tattoo removal handle, 808 diode laser
Power supply	2000w
SHR frequency	1-10HZ
808 Spot size	12*20mm, 12*35mm
Nd yag laser handle	import laser bar, with total 5 probes
Energy	10-50 J/cm2
Pulse Duration	1-10ms
Skin cooling	-21
Cooling system	Air and Close Water Circutation Treatment head with semiconductor cooling temperature adjustable
Display	12 inch TFT true color LCD touch screen
Voltage	220v/50hz or 110v/60hz
Packing size	116*60*61cm
Gross weight	90kg



15 Years Laser Beauty Device Manufacturer

Make your beauty business easier

1. Weifang Sunshinekj electronics Co., Ltd (abbreviation: MLKJ) is a highly professional and experienced manufacturer of beauty equipment in China, has 15 years in beauty industry. We have professional Production R & D department, Can provide technical support □ 2 year warranty.

Premium OEM Factories

High-performance manufacturing capacity
R&D capability for customization
Professional certifications and approvals



Product Customization Design



Production and Processing



Program Debugging



Machine Debugging



Batch Production



Fine and Safe Packaging

Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service