



Free to be yourself

A minute of painless freezing hair removal,
Feeling the comfort and coolness that
has never been seen before.

14-inch screen

2 year warranty + work 24 hours

One minute cooling -21°C



Best **TEC Cooling** system
hair remover laser



- * USA or Germany Dilas imported laser bar. (Newest and Best Quality Golden Laser Bar)
- * More than 40 million shots , more than 3 years don't need to change new handle.
- * 755 808 1064nm is for all skin and hair style hair removal

1000w , 1200w , 1600w big power for free choice
12*20 or 12*35 mm spot size

350g lightweight **handle**, more ergonomic


Available in **black or white**




High-quality
laser bars
imported from
USA

1000/1200/1600ww high power



Fast hair removal



▶ 10bar
1000W




▶ 12bar
1200W



▶ 16bar
1600W

For all skin types




755nm / 808nm / 1064 nm




Cooling of 755 808 1064 diode laser hair removal machine

- * turn on machine, immediately cooling , immediately hair removal
- * good cooling longer handle lifespan
- * real ice painless hair removal effect

Treatment effect



12*20mm or 12*35mm spot size suitable for all parts hair removal
3---5 sessions reach ideal hair removal result , permanent hair removal



- 14 inch big touch screen of 755 808 1064 diode laser hair removal machine
- * man and woman choice for all parts hair removal
 - *Any language you wanted
 - *Free logo service
 - *customized screen service