



1. Three wavelength diode laser.
2. Golden laser bars from German or USA improved 40,000,000+ longer lifespan.
3. Safe, Fast, Painless, permanent
4. For all skin and hair style hair removal
5. 14.0 inch true color touch screen with 6 skin colors.
6. Two models can be selected, single handle or double handle



Most advanced TEC cooling system inside this machine to deliver a super good cooling results.
24 hours working time. Sapphire material. Actual Pain-free experience!

One minute cooling -21°C

Best TEC Cooling system
hair remover laser

- 1) 10 bars/ 12 bars/ 16 bars, 10-40 million shots.
- 2) Big spot 12*24mm / 12*36mm for fast hair removal.
- 3) Frozen ice cooling sapphire crystal, for painless hair removal.
- 4) High quality gold laser bar imported from Germany.
- 5) Medical handle as low as 350g easy to treatment.

1 minute
and you're done

35mm {

FASTER

Screen system

1. 14 inch big touch screen, humanized designed.
2. Male/female, body parts, 6 skin types for option.
3. 6 language systems, native language + whatever you want.
4. Free OEM program service (logo/system...)

Product name	Lightsheer Diode Laser 755+808+1064nm
Laser bar	Germany/USA imported, 8 bar, 10 bar, 12 bars, 16 bars
Wavelength	755nm 808nm 1064nm/ or single 808nm
Energy	0-150J/cm ²
Frequency	1-10hz
Pulse width	40-800ms
Operative LCD Interface	14.0 TFT True Color LCD touch screen

Cooling system	Air and Close Water Circutation Treatment head with semiconductor cooling temperature adjustable
Power	4000W
Electrical source	220V/50HZ or 110V/60HZ
N.W/ G.W	62kg/78kg

About mingliangkj:

WeiFang Mingliang is a professional beauty device company specializes in the development,Our company has 14 years in beauty industry ,and pay attation to the quality and price of products .production, sales and service of beauty equipment, including IPL series, 808nm diode laser series, multi-function beauty machine series, Elight series, Q Switch ND YAG laser series, Cavitation vacuum RF series, Single RF series and so on.we are especially focus on portable and vertical multifunctional beauty machine, are very popular in the beauty machine market.

Company Strength

The collage consists of four images arranged in a 2x2 grid, each with an icon and text overlay:

- Top Left:** A person working on a computer with a circuit board on the desk. Icon: Two people. Text: **50+** R&D team.
- Top Right:** Two people working together on a project. Icon: Handshake. Text: **500+** Cooperation Case.
- Bottom Left:** A factory setting with workers and machinery. Icon: Calendar. Text: **17 Days** Product deadline.
- Bottom Right:** A row of beauty machines on wheels. Icon: Globe. Text: **100+** Export country.


- >>OEM your own machine body.
- >>OEM your own program software.
- >>OEM pretty machine colors.
- >>Print your LOGO on machine body and screen.

Custom Section


CASE design



LOGO design



SYSTEM design




Certificate Display

CERTIFICATE DISPLAY



Packing

safe and thoughtful packaging service



Built-in prevention Collision bubble

Two-year warranty for machine handle as consumable,
one-year warranty.

OEM&ODM service for distributor.

Tight packaging

Training: (Video+Manual+Online Service)
Can Keep You Operate It Easily.



3-5 business days delivery

Professional guide:
24 Hours Professional online customer service,

