



- 1.3 waves 1064nm +Diode 808nm + 755nm, 3 in 1 all skin type can be used  
High power: NEW 1000W/ 1200W /1600W handle for options
- 2.Big spot size 12\*20mm, 12\*35mmface & body hair removal solution!
- 3.Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
- 4.Long life: double life support guarantee

**Handle of diode laser hair removal machine:**





**Different spot size** Handle color optional



**Safe Use On All Skin Type**

Alex 755nm for light skin  
Speed Diode 808nm for yellow skin  
Long Pulse 1064nm for darker skin

# Treatment effect



## Software & Program

1. Male & femal, 6 body parts, 6 skin types setting smart software
2. 6 language options, your native language can be added
3. OEM program service available(logo/system...)



