

Stronger, Faster

×

14-inch HD 4K Screen



Product Description

----Safe permanent hair removal on all skin types; no pain, very comfortable during the treatment;

----Suitable for any unwanted hairs on areas like face, arms, armpits, chest, back, bikini, legs...

----You could choose 808nm or 755+808+1064nm handle, Treat all colour hair effectively.

----1 minute fast cooling. Frequency 1-10hz. Treatment fastly!!!-

NEW



NEW



NEW



For more product details, please send us inquiry now .
We have 24-hour online communication service!

Free to be yourself

1 minute of painless hair removal,
Feeling the comfort and coolness that
has never been seen before.

14-inch screen

2 year warranty + work 24 hours



About Handle

350g lightweight **handle**, more ergonomic

Available in **black or white**





35mm



before

after

Large 3.5cm spot size, and integration of three wavelengths
Hair removal is quick, efficient and convenient





One minute cooling -21°C




Best **TEC Cooling** system
hair remover laser

laser rod imported from America. 10million—40million shots. Permanent hair removal



755 808 1064nm for all different skin hair removal .

- * USA Coherent imported laser bar.
- * 10-40 million shots (Newest Golden Laser Bar)
- * 10bar 12bar 16bar ,to choose from .
- * 1000w/1200w/1600w big power,



Treatment Effect

Before and After

Treatment comparison



