



13 anos de fábrica profissional!

10 anos Alibaba Golden Fornecedor!

5.0 Alto ponto de feedback do Alibaba! Serviço ODM/OEM!



14 inch
Screen Size



3000W
Input power



1-10HZ
Frequency



755nm 808nm 1064nm
Wavelength

>>Machine 14-inch large screen

>>3000W machine power to ensure the treatment effect

>>1-10Hz fast frequency and 15x25mm spot size ensures faster, smooth and effective treatment

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>Easy to operate Data Base system,can do the treatment by man mode and women mode, também diferentes partes do corpo. Para garantir um parâmetro seguro para o novo operador

>> Long life time More than 40 millions shots

Módulo de diodo laser de ouro, anti-alta temperatura.



>>Double filters to keep water quality

>>TEC cooling Air Cooling Water CoolingSemiconductor cooling

, mantenha a temperatura da água 26 °C-28 °C o tempo todo

>>24 hours continues work

>>Easy to operate Data Base system,can do the treatment by man mode and women mode,

também diferentes partes do corpo. Para garantir um ambiente seguro parâmetro para novo operador

>> Long life time More than 40 millions shots

Módulo de diodo laser de ouro, anti-alta temperatura.



Descrição do produto



A alça pesa apenas 350g, o que é conveniente para o tratamento

>>12 pieces or 16 pieces Germany Dilas original laser bars,600/800W/1200W power is capable to meet all demands of hair removal on the whole body;

>>1-10Hz fast frequency and 15x25mm spot size ensures faster, smooth and effective treatment;

>>Super sapphire touch cooling system ensures painless and comfortable treatment.

1)Tela pode adicionar seu próprio logotipo

2) tela com 6 idiomas

3)diferentes partes do corpo para escolher, um toque pode iniciar o seu tratamento!

CERTIFICATE

Certificate – Сертификат – 證明書 – Certificat – 증명서 – 証明書

Dr. Peter Vardy
Head of Christianity Study

