

NEW ARRIVAL



C2 Vertical Magnetic Thin

- ◆ **BUILD MUSCLE**
- ◆ **BURN FAT**

The newest air-cooling system
4 working handles





Vantagens da máquina de emagrecimento

1. Deitado por 30 minutos = 5,5 horas de exercício
2. Tratamento não invasivo sem desconforto, indolor
3. 19% de Redução de Gordura
4. Aumento de 16% na massa muscular
5. Sem tempo de inatividade
6. Melhore a confiança na aparência
7. Resultados visíveis e de aparência natural
8. Operação simples e tempos de tratamento curtos



Operation



Hamstring



Thighs



Abdomen




Arm



Hip



BEFORE AND AFTER



Before

5 Weeks After

8 Weeks After



Before

After



Before

After





 **Smayer**

15 Years Laser Beauty Device Manufacturer

Make your beauty business easier



Projeto de P&D Equipes de projeto e desenvolvimento dos Estados Unidos e Alemanha.

Linha de produção madura certificada e testada e equipe de inspeção de qualidade e linha de produção padrão.

O que mais nos preocupa é que os produtos de alta qualidade sejam o padrão para manter a cooperação.

Our Factory Show

**15 Years
Professional Factory**






**Professional
Production R & D
Sales Department**



**Mass production
of various styles**



**Real Feedback
of Alibaba**



Certificate display

Certificate is the best guarantee for you



Embalagem e envio

Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service

Tempo de entrega: 3-7 dias úteis após o pagamento

Forma de entrega: por via aérea, DHL, TNT, Fedex ou agente de entrega

Antes da entrega: teste de vídeo e fotos

Garantia após a entrega: 2 anos, oferece serviço pós-venda impecável

CONTACT US



In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV SFDA ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month

Perguntas frequentes

Q1. Quando a máquina será entregue para mim?

R: Faremos a entrega em até 7 dias após o pagamento.

Q2. Você fornece serviços de treinamento?

R: Fornecemos serviços de treinamento online, após o treinamento, enviaremos um certificado.

Q3. E se a máquina quebrar?

R: Sob garantia, enviaremos componentes gratuitos, a máquina é projetada por módulo.

Cada peça é fácil de substituir.

Q4. Como fazer o pagamento?

R: aceitamos pagamento de garantia comercial alibaba, transferência bancária, western union, visa.

Q5. E a garantia?

A: 2 anos para o corpo da máquina, 1 ano para o punho. Sob garantia, se a alça quebrar, enviaremos uma nova alça.