



**MULTI-FUNCTION
4 IN 1 MACHINE**

Descrição dos produtos

Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



Laser de diodo Elight / IPL Nd yag RF, 4 alças diferentes em uma máquina

1. Este 4 em 1 RF Elight IPL laser de diodo laser nd yag tem acessórios de alta qualidade, poder forte * uso do punho do laser de diodo barra de laser dourada importada da Alemanha, 50 milhões de tiros de tempo.


* O cabo Elight usa lâmpada do Reino Unido, 150.000-30.000 tiros

*Nd yag:Equipe 5 cabeças de tratamento, 1064nm/532nm ajustável, 1064nm/532nm fixo, 1320nm fixo.


*RF: Equipa 6 cabeças de tratamento, usadas para remoção de rugas, lifting de pele, aperto de pele.

Three wavelengths work together

12*24mm² {




12*35mm² {



**Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.**





ELIGHT HANDLE

specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15*50mm, 15*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm
Standard: 420/510/680nm



Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

HANDLE-ADVANTAGE



specification

- Multiple choice
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction
1320nm for skin rejuvenation
- Indicator light
Dot position is the indicator light position dot more accurate



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm



Fixed
1064nm




Before and after



RF Handl e

12mm eye-rf	28mm small tripolar rf head	36mm big tripolar rf head	14mm bipolar-rf	28mm bipolar-rf	38mm bipolar-rf

BEFORE & AFTER



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

Advantage



808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow embroidery removal, freckle removal.

RF handle

Skin tightening/lifting, wrinkle removal on body

**Fou mais informações sobre introdução do produto e detalhes de uso, consulte-nos imediatamente
On-line 24 horas!!**