

Описание продуктов

Real shot display




Диодный лазерElight/IPLNd yagRФ, 4 разные ручки в одной машине

1.Этот 4 в 1 RFElight IPLдиодный лазериd уад лазер имеет высококачественные аксессуары, мощную мощность*Диодный лазерный манипулятор использует золотой лазерный стержень, импортированный из Германии, время 50 миллионов снимков.

- * Ручка Elight использует британскую лампу, 150 000-30 000 выстрелов *Nd yag: Оборудуйте 5 лечебных головок, регулируемые 1064нм/532нм, фиксированные 1064нм/532нм, фиксированные 1320нм.
- * RF: Оборудуйте 6 лечебными головками, используемыми для удаления морщин, подтяжки кожи.

Three wavelengths work together



Suitable for hair removal on any part of the body, Such as arms,legs, face, armpits, back, chest, etc.





















ELIGHT HANDL



specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15*50mm,15*35mm
- Long shots time: 150000 times
- 7 filters:420/510/560/585/640/690/755nm Standard:420/510/680nm



Before and after



HANDLE-ADUANTAGE



specification

Multiple choice

fixed 1064nm &532nm, adjustable 1064nm&532nm(adjust 1-7 mm)for tattoo removal, nevus removal etc.

- Multifunction
 1320nm for skin rejuvenation
- Indicator light

Dot position is the indicator light position dot more accurate











Black doll 1320nm

Fixed 532nm

Fixed 1064nm











Before and after





BEFORE & AFTE (1)







Skin rejuvenation



Freckle removal







Vascalar



Lose weight



F или дополнительную информацию о представлении продукта и деталях использования, пожалуйста, немедленно свяжитесь с нами Онлайн 24 часа!!