



Описание продуктов

# Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



## **Диодный лазер Elight/IPL Nd yag RF, 4 разные ручки в одной машине**

1. Этот 4 в 1 RF Elight IPL диодный лазер nd yag лазер имеет высококачественные аксессуары, мощную мощность\* Диодный лазерный манипулятор использует золотой лазерный стержень, импортный из Германии, время 50 миллионов снимков.


\* Ручка Elight использует британскую лампу, 150 000-30 000 выстрелов

\* Nd yag: Оборудуйте 5 лечебных головок, регулируемые 1064нм/532нм, фиксированные 1064нм/532нм, фиксированные 1320нм.


\* RF: Оборудуйте 6 лечебными головками, используемыми для удаления морщин, подтяжки кожи, подтяжки кожи.

## *Three wavelengths work together*

**12\*24mm<sup>2</sup>** {




**12\*35mm<sup>2</sup>** {



**Suitable for hair removal on any part of the body,  
Such as arms, legs, face, armpits, back, chest, etc.**





# ELIGHT HANDLE

## specification



- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15\*50mm, 15\*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm  
Standard: 420/510/680nm



## Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

# HANDLE-ADVANTAGE



## specification

- Multiple choice  
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction  
1320nm for skin rejuvenation
- Indicator light  
Dot position is the indicator light position dot more accurate



Adjustable  
532nm



Adjustable  
1064nm



Black doll  
1320nm



Fixed  
532nm



Fixed  
1064nm



# Before and after



12mm eye-rf

28mm small tripolar rf head

36mm big tripolar rf head

14mm bipolar-rf

28mm bipolar-rf

38mm bipolar-rf

# BEFORE & AFTER **R**



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

## Advantage **e**



### 808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

### Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

### ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow, embroidery removal, freckle removal.

### RF handle

Skin tightening/lifting, wrinkle removal on body

**Для или дополнительную информацию о представлении продукта и деталях использования, пожалуйста, немедленно свяжитесь с нами  
Онлайн 24 часа!!**