



MULTI-FUNCTION 4 IN 1 MACHINE

Описание продуктов

Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



Диодный лазер Elight/IPL Nd yag RF, 4 разные ручки в одной машине

1. Этот 4 в 1 RF Elight IPL диодный лазер nd yag лазер имеет высококачественные аксессуары, мощную мощность* Диодный лазерный манипулятор использует золотой лазерный стержень, импортный из Германии, время 50 миллионов снимков.


* Ручка Elight использует британскую лампу, 150 000-30 000 выстрелов

* Nd yag: Оборудуйте 5 лечебных головок, регулируемые 1064нм/532нм, фиксированные 1064нм/532нм, фиксированные 1320нм.


* RF: Оборудуйте 6 лечебными головками, используемыми для удаления морщин, подтяжки кожи, подтяжки кожи.

Three wavelengths work together

12*24mm² {




12*35mm² {



Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.





ELIGHT HANDLE

specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15*50mm, 15*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm
Standard: 420/510/680nm



Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

HANDLE-ADVANTAGE



specification

- Multiple choice
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction
1320nm for skin rejuvenation
- Indicator light
Dot position is the indicator light position dot more accurate



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm



Fixed
1064nm



Before and after



RF Handl e

12mm eye-rf	28mm small tripolar rf head	36mm big tripolar rf head	14mm bipolar-rf	28mm bipolar-rf	38mm bipolar-rf

BEFORE & AFTER



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

Advantage



808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow embroidery removal, freckle removal.

RF handle

Skin tightening/lifting, wrinkle removal on body

**Для дополнительной информации о представлении продукта и деталях использования, пожалуйста, немедленно свяжитесь с нами
Онлайн 24 часа!!**